MyTaq™ HS Red DNA Polymerase

Shipping: On dry / blue ice Catalog numbers:

BIO-21114: 250 Units

Batch No.: See vial BIO-21115 : 1000 Units Concentration: $5U/\mu L$ BIO-21116 : 2500 Units

Store at -20°C



A Meridian Life Science® Company

Storage and stability:

The MyTaq HS Red is shipped on dry/blue ice. On arrival store at -20 °C for optimum stability. Repeated freeze/thaw cycles should be avoided.

Expiry:

When stored under the recommended conditions and handled correctly, full activity of the kit is retained until the expiry date on the outer box label.

Safety precautions:

Please refer to the material safety data sheet for further information.

Unit definition:

One unit is defined as the amount of enzyme that incorporates 10 nmoles of dNTPs into acid-insoluble form in 30 minutes at 72 $^{\circ}\text{C}.$

Notes:

Research use only.

Description

MyTaq™ HS Red DNA Polymerase is a high performance PCR product that is powered by antibody-mediated hot-start, specifically designed for fast, highly-specific, hot-start PCR. MyTaq HS Red does not possess polymerase activity during the reaction set-up, thus reducing non-specific amplification including primer-dimer formation. The advanced formulation of MyTaq HS Red allows fast cycling conditions, considerably reducing the reaction time without compromising PCR specificity or yield. This new enzyme from Bioline is supplied with 5x MyTaq buffer system, a proprietary formulation that saves time and delivers superior results, as it contains dNTPs, MgCl₂ and enhancers at optimal concentrations which removes the need for optimization.

The specially designed MyTaq Red formulation does not interfere with the PCR and enables users to load samples directly onto a gel after the PCR without the need to add loading buffer.

MyTaq HS Red only requires the addition of template, primers and water, reducing the risk of pipetting errors and contamination as well as reducing the set-up time.

Components

	250 Units	1000 Units	2500 Units
MyTaq HS Red DNA Pol	1 x 50 μL	1 x 200 μL	2 x 250 μL
5x MyTaq Red Reaction Buffer	2 x 1 mL	8 x 1 mL	14 x 1.5 mL

Standard MyTaq HS Red Protocol

The following protocol is for a standard $50\mu l$ reaction and can be used as a starting point for reaction optimization. Please refer to the Important Considerations and PCR Optimization section.

PCR reaction set-up:

5x MyTaq Red Reaction Buffer	10 μL
Template	as required
Primers 20 μM each	1 μL
MyTaq HS Red DNA Polymerase	1 μL
Water (ddH₂O)	up to 50 μL

Standard PCR cycling conditions:

Step	Temperature	Time	Cycles
Initial denaturation	95 °C	1 min	1
Denaturation	95 °C	15 s	
Annealing*	User determined	15 s	25-35
Extension*	72 °C	10 s	

^{*} These parameters may require optimization, please refer to the PCR optimization section if needed

Colony PCR Protocol

MyTaq HS Red can be used for amplification of plasmid DNA directly from liquid cultures or from colonies on agar plates:

- From liquid culture: up to 8 μL of the overnight culture can be directly added to the final reaction mix.
- From colonies: we recommend using a sterile tip to stab the colony and resuspend it directly in the 50 μ L reaction mix.

Recommended cycling conditions for colony PCR of fragment up to 1 kb.

Step	Temperature	Time	Cycles
Initial denaturation	95 °C	1 min	1
Denaturation	95 °C	15 s	
Annealing*	User determined	15s	25-35
Extension*	72 °C	10 s	

^{*} These parameters may require optimization, please refer to the PCR optimization section if needed.

Multiplex PCR Protocol

MyTaq HS Red is suitable for multiplex PCR but adjustment of the annealing temperature and extension time may be required.

- Annealing temperature: We suggest using 55 °C as a starting annealing temperature. If further optimization is required we recommend using a temperature gradient to determine the optimal annealing temperature needed for the multiplex PCR.
- Extension time: since multiplex PCR require generally longer extension step, we suggest starting with a minimum of 90 s and increasing it if required.

Recommended cycling conditions for multiplex PCR of fragment up to 1 kb.

Step	Temperature	Time	Cycles
Initial denaturation	95 °C	1 min	1
Denaturation	95 °C	15 s	
Annealing*	User determined	15 s	25-35
Extension*	72 °C	90 s	

^{*} These parameters may require optimization, please refer to the PCR optimization section if needed.

Important Considerations and PCR Optimization

The optimal conditions will vary from reaction to reaction and are dependent on the template/primers used.

5x MyTaq Reaction Buffer: The 5x Reaction Buffer comprises of 5 mM dNTPs, 15 mM MgCl₂, stabilizers and enhancers. The concentration of each component has been extensively optimized, reducing the need for further optimization. Additional PCR enhancers such as Betaine etc. are not recommended.

Primers: Forward and reverse primers are generally used at the final concentration of 0.2-0.6 μ M each. As a starting point, we recommend using 0.4 μ M final concentration (*i.e.* 20 pmol of each primer per 50 μ L reaction volume). Too high a primer concentration can reduce the specificity of priming, resulting in non-specific products.

When designing primers we recommend using primer-design software such as Primer3 (http://frodo.wi.mit.edu/primer3) or visual OMPTM (http://dnasoftware.com) with monovalent and divalent cations concentrations of 10 mM and 3 mM respectively. Primers should have a melting temperature (Tm) of approximately 60 °C

Template: The amount of template in the reaction depends mainly on the type of DNA used. For templates with low structural complexity, such as plasmid DNA, we recommend using 50 pg-10 ng DNA per 50 μ L reaction volume. For eukaryotic genomic DNA, we recommend a starting amount of 200 ng DNA per 50 μ L reaction, this can be varied between 5 ng-500 ng. It is important to avoid, using template resuspended in EDTA-containing solutions (e.g. TE buffer) since EDTA chelates free Mg²⁺.

Initial Denaturation: The initial denaturation step is required to activate the enzyme and fully melt the template. We recommend 1 minute of initial denaturation at 95 °C, however for more complex templates such as eukaryotic genomic DNA, longer initial denaturation times of up to 3 minutes might be required.

Denaturation: Our protocol recommends a 15 s cycling denaturation step at 95 °C which is also suited to GC-rich templates. For low GC content (40-45%) amplicons, the denaturation step can be decreased down to 5 s, however for templates containing high GC content of >55%, a longer denaturation time of 15 s is recommended.

Annealing temperature and time: The optimal annealing temperature is dependent upon the primer sequences and is usually 2-5 °C below the lower Tm of the pair. We recommend starting with a 55 °C annealing temperature and, if necessary, to run a temperature gradient to determine the optimal annealing temperature. Depending on the reaction the annealing time can also be reduced down to 5 s.

Extension temperature and time: The extension step should be performed at 72 $^{\circ}$ C. The extension time is dependent on the length of the amplicon and the complexity of the template. An extension time of 10 s is recommended for amplicons up 1 kb or up to 5 kb in the case of a low complexity template such as plasmid DNA. For amplification of fragments over 1 kb from high complexity template, such as eukaryotic genomic DNA, longer extension times are recommended. In order to determine the fastest cycling protocol, we suggest increasing the extension time up to 30 s/kb.

Troubleshooting Guide

Problem	Possible Cause	Recommendation	
	Missing component	- Check reaction set-up and volumes used	
	Defective component	Check the aspect and the concentrations of all components as well as the storage conditions. If necessary test each component individually in controlled reactions	
No PCR	Enzyme concentration too low	- Increase enzyme quantity to up to 5U/50 μL reaction	
product	Cycling conditions not optimal	 Decrease the annealing temperature Run a temperature gradient to determine the optimal annealing temperature Increase the extension time, especially if amplifying long target Increase the number of cycles 	
	Difficult template	- Increase the denaturation time	
	Excessive cycling	- Decrease the number of cycles	
Smearing	Extension time too long	- Decrease the extension time	
or	Annealing temperature too low	- Increase the annealing temperature	
Non-Specific	Primer concentration too high	- Decrease primer concentration	
products	Contamination	 Replace each components in order to find the possible source of contamination Set-up the PCR and analyze the PCR product in separated areas. 	

Technical Support

If the troubleshooting guide does not solve the difficulty you are experiencing, please contact your local distributor or our Technical Support with details of reaction set-up, cycling conditions and relevant data.

Email: tech@bioline.com

Associated Products

Product Name	Pack Size	Cat. No.
Agarose	500 g	BIO-41025
Agarose tablets	300 g	BIO-41027
HyperLadder™ 1kb	200 Lanes	BIO-33025

TRADEMARK AND LICENSING INFORMATION

- 1. Notice to Purchaser: Licensed under U.S. patent numbers 5,338,671 and 5,587,287 and corresponding patents in other countries
- 2. HyperLadder and MyTaq are Trademarks of Bioline Reagents Ltd

Bioline GmbH GERMANY Bioline Reagents Ltd UNITED KINGDOM Bioline (Aust) Pty. Ltd AUSTRALIA Bioline USA Inc. Bioline France Meridian Bioscience Asia Pte Ltd **FRANCE** SINGAPORE Tel: +44 (0)20 8830 5300 Tel: +1 508 880 8990 Tel: +49 (0)337 168 1229 Tel: +61 (0)2 9209 4180 Tel: +33 (0)1 42 56 04 40 Tel: +65 6774 7196 Fax: +44 (0)20 8452 2822 Fax: +1 508 880 8993 Fax: +49 (0)3371 68 1244 Fax: +61 (0)2 9209 4763 Fax: +33 (0)9 70 06 62 10 Fax: +65 6774 6441