# HyperLadder™ 100bp

Shipping: Room Temperature

Catalog numbers

Batch No.: See vial BIO-33056: 100 Lanes

BIO-33029: 200 Lanes

BIO-33030: 500 Lanes



A Meridian Life Science® Company

Storage and stability:

HyperLadder 100bp is shipped at room temperature. On arrival store at -20°C for optimum stability.

Expiry:

When stored under the recommended conditions and handled correctly, band quality is retained until the expiry date on the outer box label.

#### Safety precautions:

Please refer to the material safety data sheet for further information.

#### Quality control specifications:

HyperLadder 100bp is extensively tested quality prior to release.

#### Notes:

Research use only.

#### **Features**

- Easy size determination
- 10 sharp and bright bands from 100bp 1013bp
- Ready-to-use format
- Optional mass determination

# **Description**

HyperLadder 100bp (formerly HyperLadder IV) is a ready-to-use molecular weight marker, especially designed for easy size determination. This ready-to-use format reduces handling steps and saves time; simply transfer HyperLadder IV from the vial to the gel.

HyperLadder IV produces a pattern of 10 regularly spaced bands, ranging from 100bp to 1013bp. To allow easy identification and orientation, the 1000/1013bp band has the highest intensity.

Optionally the HyperLadder can be used for mass determination.

A 5x sample-loading buffer is supplied for your convenience. Under no circumstances should it be used to dilute/load ladder.

## Components

Product	100 Lanes	200 Lanes	500 Lanes
HyperLadder <sup>™</sup> 100bp	500µl	2 x 500µl	5 x 500µl
5x Sample Loading Buffer	1ml	1ml	1ml

# **Associated Products**

Product	Pack size	Cat. No.
Agarose	100g	BIO-41026
Agarose Tablets	150g	BIO-41028


# **Protocol**

Vortex the tube and use  $5\mu l$  per lane. There is no need to add loading buffer to the ladder.

Use the 5x loading buffer for sample loading, under no circumstances should it be used for diluting the ladder. The blue dye in the loading buffer runs at approximately 4bp on a 1% gel.

### Composition

This marker is composed of a restriction digest plus one or more PCR products.



#### Notes

- 1. This product insert is a declaration of analysis at the time of manufacture
- 2. HyperLadder is a trademark of Bioline
- 3. No licence is conveyed in respect of US 6,680,378
- 4. Research Use Only

Bioline Reagents Ltd UNITED KINGDOM

Tel: +44(0)20 8830 5300 Fax: +44 (0)20 8452 2822 Bioline USA Inc. USA

Tel: +1 508 880 8990 Fax: +1 508 880 8993 Bioline GmbH GERMANY

Tel: +49 (0)33 7168 1229 Fax: +49 (0)33 7168 1244 Bioline (Aust) Pty. Ltd AUSTRALIA

Tel: +61 (0)2 9209 4180 Fax: +61 (0)2 9209 4763