# **Tetro Reverse Transcriptase**

Shipping: On Dry/Blue Ice Catalog numbers

> BIO-65050 10,000 units

> > 4 x 10.000 units

Batch No.: See vial Concentration: 200u/ul



Store at -20°C

BIO-65051

Tetro Reverse Transcriptase is shipped on dry/blue ice. All kit components should be stored at -20°C upon receipt. Excessive freeze/thawing is not recommended.

When stored under the recommended conditions and handled correctly, full activity of the kit is retained until the expiry date on the outer box label.

#### **Quality Control:**

Storage and stability:

The Tetro Reverse Transcriptase and its components are extensively tested for activity, processivity, efficiency, sensitivity, absence of nuclease contamination and absence of nucleic acid

#### Safety Precautions:

Please refer to the material safety data sheet for further information.

#### Notes:

For research use only.

### Description

Tetro Reverse Transcriptase is a highly sensitive, high stability M-MLV reverse transcriptase. Tetro Reverse Transcriptase is optimized for reverse transcription reactions using a wide range of total RNA amounts (100pg-2µg), such that long and low abundance cDNAs can be detected by amplification after cDNA synthesis.

Tetro Reverse Transcriptase is suitable for first-strand cDNA synthesis, cDNA library construction, and the production of templates for RT-PCR analysis of gene expression. Tetro Reverse Transcriptase can be used with total RNA, mRNA, in-vitro transcribed RNA or viral RNA.

#### Kit components

| Reagent                     | 10,000 units | 40,000 units |
|-----------------------------|--------------|--------------|
| Tetro Reverse Transcriptase | 50µI         | 4 x 50µl     |
| 5x Reaction Buffer          | 1.2ml        | 4 x 1.2ml    |

#### **Reaction Recommendations and Optimization**

#### **Template Quality**

- Intact, high-quality RNA is essential for the reverse-transcription
- All reagents for use with RNA must be prepared using DEPC-treated water (BIO-38030).
- The inclusion of an RNase Inhibitor can reduce template degradation and increase yield of PCR product (BIO-65027).
- Low-copy-number genes may require an increase in starting material.
- It is necessary to use a suitable RNA extraction reagent e.g., TRIsure<sup>TM</sup> (BIO-38032) or RNA Isolation Kit (BIO-52072).

#### **Primer Design and Concentration**

There are three methods for priming cDNA synthesis:

#### Oligo dT Primers

Oligo dT priming (BIO-38029) uses the poly-A tail found on the 3' end of most eukaryotic mRNAs. This ensures that the 3' end of mRNAs are represented, although long mRNAs can have their 5' ends underrepresented in the subsequent cDNA pool. Use at 0.5µM final concentration.

#### **Random Hexamers**

Random priming (BIO-38028) gives random coverage to all regions of the RNA to generate a cDNA pool containing various lengths of cDNA. Random priming is unable to distinguish between mRNA and other RNA species present in the reaction. Use at 2.0µM final concentration.

# Gene Specific Primers (GSP)

Gene specific primers are designed to generate cDNA for a specific gene of interest. It is a widely used method for performing One-Step RT-PCR when only 1 gene is under investigation. It can be useful when RNA concentrations are low. Use at 0.4µM final concentration.

A combination of Oligo dT and Random Hexamers primers can improve the reverse transcription efficiency of some mRNA templates.

# **Extension Temperature**

- Efficient reverse-transcription can be achieved at temperatures of 37°C to 45°C for 30-60 min.
- The use of higher incubation temperatures up to 48°C may increase the yield of cDNA synthesized in cases of complex RNA secondary structure. However, the yield of the majority of RNA molecules will be reduced.

#### **Tetro Reverse Transcription Protocol**

- 1. Vortex solutions and centrifuge briefly before use.
- 2. Prepare the priming premix on ice in an RNase-free reaction tube:

| Total RNA (up to 5µg) or mRNA (up to 0.5µg) |   | <i>n</i> µl |
|---|---|-------------|
| Primer*:                                    | Oligo (dT) <sub>18</sub> (10μM) <i>or</i><br>Random Hexamer (40μM) <i>or</i><br>GSP (8μM) | 1µl         |
| 10mM dNTP mix*                              |   | 1µl         |
| 5x RT Buffer                                |   | 4µl         |
| RiboSafe RNase Inhibitor                    |   | 1µl         |
| Tetro Reverse Transcriptase (200u/µl)       |   | 1µl         |
| DEPC-treated water*                         |   | to 20µl     |

<sup>\*</sup> Available separately (see Associated Products)

- 3. Mix gently by pipetting.
- 4. Incubate samples at 45°C for 30 min. If using random hexamers, incubate 10 min at 25°C followed by 45°C for 30 min.
- 5. Terminate reaction by incubating at 85°C for 5 min, chill on ice.
- 6. Store reaction at -20°C for long term storage, or proceed to PCR immediately.

This protocol is intended for use as a guide only; conditions will vary from reaction to reaction and may need optimization.

Website: www.bioline.com/ email: info@bioline.com PI-50258 V4

## **Troubleshooting**

| Problem                     | Possible Cause                                | Recommendation  |  |
|-----------------------------|---|---|--|
| No cDNA synthesis           | RNA degraded                                  | Analyze RNA on a denaturing gel to verify integrity. Ensure that all reagents are RNase-free. Use Ribosafe RNase inhibitor in the first-strand reaction (BIO-65027).  |  |
|                             | RNA contained an RT inhibitor                 | The presence of inhibitors can be determined by mixing a control RNA with some of the sample and comparing the yield with that of the original amplification. Remove inhibitors such as SDS, EDTA, formamide and pyrophosphate, by ethanol precipitation of RNA, including a 70% ethanol wash step. |  |
|                             | Reaction temperature not optimal              | Perform a temperature-gradient experiment ranging from 37-48°C.   |  |
|                             | Not enough starting RNA                       | Increase the amount of starting RNA, this can be an important factor when amplifying low-copy genes from total RNA.   |  |
|                             | RNA had high secondary structure              | Prior to reaction set-up, denature RNA with primers. Raise the temperature of the RT step, up to a maximum of 48°C (for short amplicons).   |  |
|                             | Insufficient product                          | Increase reverse transcription step to 60 minutes   |  |
| Poor Specificity in PCR     | Non-specific annealing of primers to template | Use gene-specific primers rather than Oligo dT or random hexamers in RT reaction. Increase the annealing temperature in PCR. Check for presence of pseudogenes. Set up reactions on ice.  |  |
|                             | Primer dimers                                 | Redesign primers to prevent self-annealing.   |  |
|                             | Genomic DNA contamination                     | Treat RNA with DNase I and re-purify. If possible, use intron-spanning primers in PCR.  |  |
| Product in no-RTase control | Template contaminated with DNA                | Treat samples with DNase I.   |  |

# **Technical Support:**

If the troubleshooting guide does not solve the difficulty you are experiencing, please contact Technical Support with details of reaction setup, cycling conditions and relevant data.

Email: tech@bioline.com

# **Product Citations:**

- 1. To, K.W., et al. Mol. Can. Res. 9, 516-527 (2011).
- 2. Comerford, I., et al. Blood 116(20), 4130-4140 (2010).
- 3. Corripio-Miyar, Y., et al. Mol. Immunol. 46(10), 2098-2106 (2009).
- 4. Chen, Y., et al. Blood 114(1), 40-48 (2009).
- 5. Le, H. K., et al. Can. Immunol. Immunother. **58(10)**, 1565-1576 (2009).

# **Associated products:**

| Product Name               | Cat. No.  |
|----------------------------|-----------|
| RiboSafe RNase Inhibitor   | BIO-65027 |
| TRIsure™                   | BIO-38032 |
| ISOLATE II RNA Mini Kit    | BIO-52072 |
| Random Hexamer Primer      | BIO-38028 |
| Oligo (dT)18 Primer        | BIO-38029 |
| dNTP Mix (10mM)            | BIO-39053 |
| Tetro cDNA Synthesis Kit   | BIO-65043 |
| SensiFAST™ SYBR No-ROX Kit | BIO-98002 |
| Agarose, Molecular Grade   | BIO-41026 |

Bioline Reagents Ltd UNITED KINGDOM Bioline USA Inc.

Bioline GmbH GERMANY

Bioline (Aust) Pty. Ltd AUSTRALIA Bioline France FRANCE Meridian Bioscience Asia Pte Ltd SINGAPORE

Tel: +44 (0)20 8830 5300 Fax: +44 (0)20 8452 2822

Tel: +1 508 880 8990 Fax: +1 508 880 8993 Tel: +49 (0)337 168 1229 Fax: +49 (0)3371 68 1244 Tel: +61 (0)2 9209 4180 Fax: +61 (0)2 9209 4763 Tel: +33 (0)1 42 56 04 40 Fax: +33 (0)9 70 06 62 10 Tel: +65 6774 7196 Fax: +65 6774 6441