


## One Click® Water Determination

### Simple and Secure



#### One Click® User Interface

The large color touchscreen features the simple-to-understand and fast-to-operate One Click user interface. All users have their own homescreen with dedicated shortcut buttons to directly start specific routine tasks.




#### Solvent Manager

Users are protected against contact with chemicals using the Solvent Manager liquid handling device. It is fully controlled by the titrator and allows reagents to be removed with minimal user interaction.



#### Diaphragmless cell

The user has the choice of a current generator cell with or without a diaphragm. For most samples the cell without a diaphragm offers more convenience, as no catholyte is needed. Due to the unique design, the diaphragmless cell even meets the demands of very precise low water content determinations.



#### PC software LabX® Titration

Full control over the titrator, secure database archiving, networking, efficient result management with control charts and individual search filters or rapid and clear method development are only some aspects of the powerful LabX Titration PC software.



#### C20S Compact KF Coulometer

The Routine Karl Fischer Coulometer combines fast and precise water content determination with very simple operation. Samples from any industry with water content between 1 ppm and 5% are analyzed in single measurements or as sample series with statistics.

- Color touchscreen with One Click user interface
- User interface in English, German, French, Spanish, Italian, Chinese, Russian, Polish, Korean and Portuguese
- Generator cell with or without diaphragm available
- Connectivity to LabX Titration PC software via USB or direct Ethernet (TCP/IP)
- Common USB printers or the USB-P25 compact printer for secure and GLP-compliant paper recording

| Functionality Overview |   | C20S |
|------------------------|---|------|
| Application            | Liquid and gaseous samples  | •    |
|                        | Resources and method database   | •    |
|                        | User specific Home Screen with Shortcuts  | •    |
| One Click® Titration   | Shortcuts per user  | 12   |
|                        | Reagent Control   | –    |
|                        | Solvent Manager   | •    |
| Hot Plug & Play        | USB printer DIN-A4/US letter / USB-P25 compact printer  | •    |
|                        | USB memory stick / USB barcode reader / USB hub   | •    |
|                        | LevelSens   | –    |
| Security               | LogStraight Fingerprint Reader  | –    |
| Automation             | Stromboli oven sample changer   | –    |
| User Management        | Expert – Routine  | –    |
|                        | Preprogrammed METTLER TOLEDO methods  | –    |
|                        | KF standard methods   | 1    |
| Methods and series     | Max. number of user methods   | 5    |
|                        | External extraction/dissolution dedicated method type   | –    |
|                        | Bromine Index determination   | –    |
|                        | Number of method functions  | 12   |
|                        | Number of samples per series  | 120  |
|                        | Autostart: automatic start after sample addition  | •    |
|                        | Max. number of calculations/results per method  | 3    |
|                        | User defined calculations   | –    |
|                        | Statistics of up to 120 samples including result limits check, outlier test   | •    |
|                        | Extra statistical functions as max. srel and multiple determinations  | –    |
|                        | Recalculation of results  | •    |
|                        | Data in Online Screen: time, water count, potential, stirrer speed, consumption, drift                                  | •    |
|                        | Data in measured value table: time, water count, potential, consumption, drift  | •    |
|                        | Curves in Online Screen and printout: E-C, E-t, C-t, drift-t, water count-t   | •    |
| Task list              | Number of tasks queued  | –    |
| Online Help            | Context sensitive parameter explanation   | •    |
| Data Export            | USB, Ethernet, PDF  | •    |
| Languages              | User interface and printouts in English, German, French, Spanish, Italian, Chinese, Russian, Polish, Korean, Portuguese | •    |
| PC software            | LabX® Titration Express and Server  | •    |
|                        | Ethernet (TCP/IP)   | 1    |
|                        | COM (serial balance connection)   | 1    |
|                        | USB host (printer, memory stick, barcode reader, hub)   | 2    |
|                        | USB client (PC)   | 1    |
|                        | Integrated Magnetic Stirrer   | 1    |
|                        | Pump / Stirrer  | 1    |

- exclusive C30S functionality

## Technical Specifications

|  |   |                                     |
|--|---|-------------------------------------|
| Polarized sensor input                                   | Measurement range                                 | -2000...+2000 mV                    |
|  | Resolution / Limits of error                      | 0.1 mV / 2 mV                       |
|  | Polarization current range                        | 0...24 µA                           |
|  | Polarization current resolution / Limits of error | 0.1 µA / 1 µA                       |
| Environment  | Ambient temperature                               | +5...+40 °C                         |
|  | Atmospheric humidity                              | max. 80% at 31 °C                   |
| Titrator dimensions                                      | Width x depth x height / Weight                   | 210 x 340 x 312 mm / 3.3 kg         |
|  | Glass beaker                                      | 200 mL                              |
| Control and display                                      | Touchscreen                                       | VGA 5.7" color TFT, 640 x 480 pixel |
| Application<br>(depends on the environment and handling) | Drift measurement                                 | Online, < 2 µg/min                  |
|  | Measurement range                                 | 10 µg...200 mg                      |
|  | Resolution  | 0.1 µg                              |
|  | Repeatability                                     | 0.3% at >1 mg H <sub>2</sub> O      |



**Mettler-Toledo GmbH, Analytical**  
CH-8603 Schwerzenbach, Switzerland  
Tel. +41 44 806 77 11  
Fax +41 44 806 72 40

Subject to technical changes  
© 03/2016 Mettler-Toledo GmbH  
Printed in Switzerland 30265579A  
Marketing Titration / Global MarCom Switzerland, 1688 AK/MD

[www.mt.com/titration](http://www.mt.com/titration)

For more information

One Click is a Registered Trademark of METTLER TOLEDO in Switzerland, the European Union, Russia and Singapore.  
LabX is a Registered Trademark of METTLER TOLEDO in Switzerland, USA, China, Germany and a further 13 countries.