WaterSep BioSeparations

BioProducer12 Steamer

Autoclavable,Gamma Irradiated Hollow Fiber Cartridges. Easy and Reliable



The BioProducer12 Steamer is an autoclavable, gamma irradiated, re-usable hollow fiber cartridge made with WaterSep's low binding, glycerin free modified PES membrane (m-PES) which provides a process flux and product recovery that exceeds other cross flow devices. BioProducer Steamer cartridges are ready to use without any pre-rinse or sanitization, or can be autoclaved as is, or in an assembly.

WaterSep's **BioProducer Steamer** cartridges are true plugand-play, crossflow filtration at its best - connect, use and discard, or clean with 0.5 - 1.0 M NaOH and re-use.


The BioProducer Steamer production scale hollow fiber cartridges are designed for biopharmaceutical concentration/diafiltration and clarification unit operations. With a membrane surface area of 1.3 m², the BioProducer12 Steamer HF cartridge can accommodate batch volumes between 25L and 100L with a typical permeate flow rate of 39-78 L/hr for ultrafiltration applications and 26-52 L/hr for microfiltration applications.

The **BioProducer12 Steamer** hollow fiber cartridge products are in most applications superior to conventional techniques, like centrifugation, depth filters or UF/DF cassettes.

The **BioProducer12 Steamer** hollow fiber cartridges are offered with the same path length as other high performance WaterSep HF cartridges, which makes scale up easy and predictable.

The **BioProducer12 Steamer** cartridges offer a true single use crossflow system that is:

- Gamma irradiated ready to use no rinse.
- Humectant free 80% less extractables, see Extractables Flush Analysis.
- No time-consuming cleaning validation.
- At a reduced risk of cross contamination.
- Consistent, membrane performance batch to batch.
- Self-contained and requires no assembly, no expensive hardware and no membrane installation.
- Significantly less expensive than traditional cassette formats.
- Every HF cartridge is integrity tested and has an individual lot number for easy traceability.



TYPICAL APPLICATIONS

- **Cell perfusion** with an ATF System or conventional pump operated crossflow perfusion system. Bioreactor and **Steamer** hollow fiber assembly can be autoclaved simultaneously.
- Small to medium scale **aseptic operations** or any other sterile application requiring autoclaving.
- Rapid clarification of small volumes of **cell culture**, **fermentation solutions and virus/vaccine suspensions**.
- Concentration and diafiltration of monoclonal antibodies, enzymes, blood components and other proteins.



BioProducer12 Steamer

ORDERING INFORMATION

Place orders online at <u>www.watersep.net</u>, or email purchase orders to <u>orders@watersep.net</u>, sales order confirmation to follow.

BioProducer12 Steamer HF Cartridge

0.5 mm ID* AU XXX 05PRO12 S3, 3/pkg AU XXX 05PRO12 S6, 6/pkg *Not available in 0.65 μm

1.0 mm ID AU XXX 10PRO12 S3, 3/pkg AU XXX 10PRO12 S6, 6/pkg

2.0 mm ID* AU XXX 20PRO12 S3, 3/pkg AU XXX 20PRO12 S6, 6/pkg *Available in 10K, 50K, 300K - 750K and 0.1 μm - 0.45 μm

Replace XXX with:

003 for 3K	100 for 100K	910 for 0.1 µm
005 for 5K	300 for 300K	920 for 0.2 µm
010 for 10K	500 for 500K	945 for 0.45 µm
030 for 30K	750 for 750K	965 for 0.65 µm

SPECIFICATIONS

Cartridge Dimensions 3.50" (8.89 cm) x 15.00" (38.10 cm)

Membrane Surface Area 13.99 ft² (1.30 m²)

Molecular Weight Cut-off 3K, 5K, 10K, 30K, 50K, 100K, 300K, 500K, 750K

Membrane Pore Size 0.1 μ m, 0.2 μ m and 0.45 μ m, 0.65 μ m

Fiber ID 0.5 mm, 1.0 mm, 2.0 mm

MATERIALS OF CONSTRUCTION

The cartridge as assembled is USP Class VI compliant.

Membrane: Modified Polyethersulfone (m-PES)


Housing: White Polysulfone

Encapsulant: Epoxy

CONNECTIONS

Feed/Retentate 1.5" TC Standard Sanitary Fitting

Permeate 1.0" TC Standard Sanitary Fitting



Technical Support: <u>experttalk@watersep.net</u> or call 508.970.0089 x 204.

WaterSep BioSeparations Corporation 450 Donald Lynch Blvd. Marlborough, MA 01752

www.watersep.net 508.970.0089 x204



ISO 9001:2015 Certified Registration #19.8073