Document id: AUDIS24PS rev 1

Discover24 Steamer

Autoclavable, Gamma Irradiated Hollow Fiber Cartridges Easy and Reliable



The Discover24 Steamer is a small scale autoclavable, gamma irradiated, re-usable hollow fiber cartridge, made with WaterSep's low binding, glycerin free modifed PES membrane. The Discover24 Steamer is ready to use without any pre-rinse or sanitization, or can be autoclaved as is, or in an assembly.

WaterSep's Steamer Line is true plug-and-play, at its best connect, use and discard, or clean with 0.5-1.0 M NaOH and re-use. The extractable level is approximately 80x less than a glycerin conditioned membrane, see Extractables Flush Analysis. After a quick buffer conditioning the HF cartridge is ready to be used or autoclaved. For autoclaving directions, please see the Steamer Validation Guide at www.WaterSep.net.


The Discover Steamer hollow fiber cartridge products are in most applications superior to conventional techniques, like centrifugation, depth filters, stirred cells or dialysis bags.

The Discover24 Steamer cartridges have a membrane surface area of 107 cm², providing a typical permeate flow rate of 32 to 64 ml/hr for ultrafiltration applications and 21 to 43 ml/hr for microfiltration applications. The recommended sample volume for the Discover24Steamer is 50-400 ml.

The fluid path length for the Discover24 Steamer cartridges is similar to all other WaterSep 24" (60 cm) hollow fiber cartridges making scale up and scale down easy and predictable.

The Discover24 Steamer cartridges offer true ready to use crossflow filtration that is:

- Autoclavable as is or in an assembly.
- Gamma irradiated ready to use no rinse.
- · Humectant free 80% less extractables. see Extractables Flush Analysis.
- Free of time-consuming cleaning validation.
- · At a reduced risk of cross contamination.
- Consistent, membrane performance batch to batch.
- · Self-contained and requires no assembly, no expensive hardware and no membrane installation.
- · Significantly less expensive than traditional cassette formats.
- Every HF cartridge is integrity tested and has an individual lot number for easy traceability.



TYPICAL APPLICATIONS

- Cell perfusion with an ATF system or conventional pump operated perfusion. (Bioreactor and Steamer hollow fiber assemblies can be autoclaved simultaneously.)
- Small aseptic/sterile UF-MF applications where autoclaving is required.
- · Rapid clarification of small volumes of cell culture, fermentation solutions and virus/vaccine suspensions.
- Concentration and diafiltration of monoclonal antibodies, enzymes, blood components and other proteins.
- Preparative biopharmaceutical work.





Discover24 Steamer

ORDERING INFORMATION

Place orders online at www.watersep.net, or email purchase orders to orders@watersep.net, sales order confirmation to follow.

Discover24 Steamer HF Cartridge

0.5 mm Fiber ID* AU XXX 05DIS24 L6, 6/pkg AU XXX 05DIS24 LD, 12/pkg *Not available in 0.65 μm

1.0 mm Fiber ID AU XXX 10DIS24 L6, 6/pkg AU XXX 10DIS24 LD, 12/pkg

2.0 mm Fiber ID* AU XXX 20DIS24 L6, 6/pkg AU XXX 20DIS24 LD, 12/pkg *Available in 10K, 50K, 300K - 750K and 0.1 μ m - 0.45 μ m

Replace XXX with:

 $\begin{array}{ccccc} 003 \; \text{for 3K} & 100 \; \text{for 100K} & 910 \; \text{for 0.1} \; \mu\text{m} \\ 005 \; \text{for 5K} & 300 \; \text{for 300K} & 920 \; \text{for 0.2} \; \mu\text{m} \\ 010 \; \text{for 10K} & 500 \; \text{for 500K}. & 945 \; \text{for 0.45} \; \mu\text{m} \\ 030 \; \text{for 30K} & 750 \; \text{for 750K} & 965 \; \text{for 0.65} \; \mu\text{m} \\ \end{array}$

050 for 50K

SPECIFICATIONS

Cartridge Dimensions 0.38" (0.95 cm) x 24.38" (61.93 cm)

Membrane Surface Area 0.115 ft² (107 cm²)

Molecular Weight Cut-off 3K, 5K, 10K, 30K, 50K, 100K, 300K, 500K, 750K

Membrane Pore Size 0.1 μ m, 0.2 μ m and 0.45 μ m, 0.65 μ m

Fiber ID 0.5 mm, 1.0 mm, 2.0 mm

MATERIALS OF CONSTRUCTION

The cartridge as assembled is USP Class VI compliant.

Membrane:

Modified Polyethersulfone (m-PES)

Housing:

White Polysulfone

Encapsulant: Epoxy


CONNECTIONS

Feed/Retentate:

Female Luer Lok Fitting

Permeate:

Female Luer Lok Fitting



Technical Support: experttalk@watersep.net or call 508.970.0089 x 204.