# perLadder™ 25b

Shipping: Room Temperature

Catalog numbers

Batch No.: See vial BIO-33057: 100 Lanes

> BIO-33031: 200 Lanes BIO-33032: 500 Lanes

HyperLadder 25bp is shipped at room temperature. On arrival store at -20°C for optimum stability.

Expiry:

When stored under the recommended conditions and handled correctly, band quality is retained until the expiry date on the outer box label.

Safety precautions:

Storage and stability:

Please refer to the material safety data sheet for further information.

Quality control specifications:

HyperLadder 25bp is extensively tested quality prior to release.

Notes:

Research use only.



**Features** 

- Easy mass determination
- 12 bands from 25bp 500bp
- Ready-to-use format
- Optional mass determination

### Description

HyperLadder 25bp (formerly HyperLadder V) is a ready-to-use molecular weight marker for size determination of DNA fragments. It is especially designed for short fragments such as apoptotic DNA oligonucleotides. This ready-to-use format reduces handling steps and saves time; simply transfer HyperLadder 25bp from the vial to the gel.

HyperLadder 25bp produces a pattern of 12 regularly spaced bands, ranging from 25bp to 500bp. To allow easy identification and orientation, the 100bp and 200bp bands have the highest intensity.

Optionally the HyperLadder can be used for mass determination.

A 5x sample-loading buffer is supplied for your convenience. Under no circumstances should it be used to dilute/load ladder.

## Components

Product	100 Lanes	200 Lanes	500 Lanes
HyperLadder <sup>™</sup> 25bp	500µI	2 x 500µl	5 x 500µl
5x Sample Loading Buffer	1ml	1ml	1ml


# **Associated Products**

Product	Pack size	Cat. No.
Agarose	100g	BIO-41026
Agarose Tablets	150g	BIO-41028

#### **Protocol**

Vortex the tube and use 5µl per lane. There is no need to add loading buffer to the ladder.

Use 5x loading buffer for sample loading, under no circumstances should it be used for diluting the ladder. The blue dye in the loading buffer runs at approximately 4bp on a 1% gel



#### Notes

- 1. This product insert is a declaration of analysis at the time of manufacture
- 2. HyperLadder is a trademark of Bioline
- 3. No licence is conveyed in respect of US 6,680,378
- 4. Research Use Only

Bioline Reagents Ltd UNITED KINGDOM

Tel: +44(0)20 8830 5300 Fax: +44 (0)20 8452 2822 Bioline USA Inc.

Tel: +1 508 880 8990 Fax: +1 508 880 8993 Bioline GmbH **GERMANY** 

Tel: +49 (0)33 7168 1229

Bioline (Aust) Pty. Ltd **AUSTRÀLIA** 

Tel: +61 (0)2 9209 4180 Fax: +61 (0)2 9209 4763