

ENZYCLEAN® II LS

Dual Enzyme Low Suds Detergent



Enzyclean® II LS Dual Enzyme Low Suds Detergent is formulated for use as a presoak or manual detergent and for use in ultrasonic washers, automatic washing equipment and evacuators. The protease and amylase enzymes in Enzyclean® II LS break down organic and inorganic soils, residues, films and deposits and ensures thorough and complete cleaning, leaving instruments moving freely and spot-free.

- **Compatible with all equipment & pump systems.**
- **Protease enzyme effectively breaks down proteins.**
- **Amylase enzyme breaks down carbohydrates.**

Dual Enzyme Low Suds Detergent

Order No.	Size	Physical Properties
Z6	1 Gallon (3.8 Liters)	Form: Liquid Color: Green Odor: Spearmint pH Concentration: 8.8

Uses	Temperature
<input checked="" type="checkbox"/> Presoak <input checked="" type="checkbox"/> Manual <input checked="" type="checkbox"/> Ultrasonic <input checked="" type="checkbox"/> Evacuators	Usage: 60°-140° F (16°-60° C) Storage: Cool, dry area

Product Features	Consumer Benefits
Dual Enzymes:	Protease enzyme effectively breaks down proteins. Amylase enzyme breaks down carbohydrates.
Highly Concentrated:	Low cost per use.
Rust Inhibitors:	Reduces likelihood of rust, tarnish or corrosion.
Rinsing Agents:	Rinses spot & residue-free.
Neutral pH:	Product safe and compatible with all metals.
Low Sudsing:	Compatible with all equipment & pump systems.
Bacteriostatic:	If stored properly, not favorable to bacterial growth.
Environmentally-Friendly:	Biodegradable, phosphate-free formula.

Procedures	Directions for Use
Directions:	Wear appropriate personal protective equipment based on safe work practices.
Presoaking, Manual and Ultrasonic Cleaning:	Add 1/8 fl. oz. (3.7 mL) to 2 fl. oz. (59 mL) per gallon (3.8 L) of warm water. Allow to soak for 2-3 minutes or longer. Discard solution after each use or when visibly soiled.
Automatic Washing Equipment:	Add 1/8 fl. oz. (3.7 mL) to 2 fl. oz. (59 mL) per gallon (3.8 L) of water. Recommended temperature range: 90°-140° F (32°-60° C).
Evacuators:	Add 1/8 fl. oz. (3.7 mL) to 2 fl. oz. (59 mL) per gallon (3.8 L) of water. Suggested maintenance twice weekly.
	Dilutions vary depending on temperature, water quality and equipment. Validated cleaning is recommended.