REF 985090

de

Test 0-90

03.23

NANOCOLOR® Sulfite 100**Methode:**

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03.23

NANOCOLOR® Sulfite 100**Method:**

Photometric determination with potassium iodate / iodide

Range:	5 – 100 mg/L SO_3^{2-}
Wavelength (HW = 5 – 12 nm):	470 nm
Reaction time:	0
Reaction temperature:	20 – 25 °C

Contents of reagent set:

- 19 test tubes Sulfite 100
- 2 test tubes with 10 mL Sulfite 100 R1
- 1 test tube with 9 mL Sulfite 100 R2
- 1 empty test tube for blank value

Hazard warning:

Reagent R2 contains sulfuric acid 30 – 51 %.

H314 Causes severe skin burns and eye damage.

P260, P280, P301+330+331, P303+361+353, P304+340, P305+351+338 Do not breathe vapors. Wear protective gloves / eye protection. IF SWALLOWED: rinse mouth. Do NOT induce vomiting. IF ON SKIN (or hair): Remove / Take off immediately all contaminated clothing. Rinse skin with water / shower. IF INHALED: Remove to fresh air and keep at rest in a position comfortable for breathing. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. For further information ask for a safety data sheet.

Preliminary tests:

If the order of magnitude of the concentration in a sample is not known, a preliminary test with QUANTOFIX® Sulfite (10 – 1000 mg/L SO_3^{2-} , REF 91306) or with VISOCOLOR® HE Sulfite SU 100 (REF 915008) rapidly gives this information. From the order of magnitude the required dilution can be calculated and prepared directly.

Interferences:

Oxidizing and reducing substances as well as starch interfere with the determination.

2 mg/L ascorbic acid simulate 1 mg/L sulfite.

The method can be used also for the analysis of sea water.

Procedure:

Requisite accessories: piston pipette with tips

Test sample	Blank value
Open test tube, add	Open empty test tube, add
200 µL (= 0.2 mL) R2, shake lightly, slowly add along wall of tube	1.0 mL R1 and
4.0 mL sample solution (<i>the pH value of the sample must be between pH 1 and 13</i>), close and mix.	200 µL (= 0.2 mL) R2, shake lightly, slowly add along wall of tube
Clean outside of test tube and measure immediately.	4.0 mL distilled water, close and mix.
	Clean outside of test tube and measure immediately.

Measurement:

For NANOCOLOR® photometers and PF-12 see manual, test 0-90.

For each analysis a blank value is required.

Sulfite contents below 5 mg/L are to be considered as negative results.

Photometers of other manufacturers:

For other photometers check whether measurement of round glass tubes is possible. Verify factor for each type of instrument by measuring standard solutions.

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MACHEREY-NAGEL GmbH & Co. KG · Valenciennes Str. 11 · 52355 Düren · Niemcy