

GemCell Plus™ U.S. Origin Human Serum AB — Xeno-Free

Catalog # 100-912

Thomas No. CHM03W512

PRODUCT DESCRIPTION:

GemCell Plus™ Xeno-Free Human Serum is collected from healthy male donors of the AB serotype at FDA-licensed facilities in the United States. This material goes through additional rigorous testing and contains therapeutic-grade human thrombin. All donor units are tested for viral markers and found to be non-reactive.

- This product was developed without the use of animal material.
- 0.1 µm sterile-filtered

Precautions: Non Hazardous Ingredients: Blood Plasma (Human)

The product contains no substances that are considered hazardous to health at their given concentration. Although the human plasma used to manufacture this product has been tested and shown to be negative for infectious agents, no test is 100% accurate. All materials derived from human plasma should be handled as if capable of transmitting infection.

INSTRUCTIONS FOR USE

Receipt and Storage:

Ships on dry ice. Store frozen at -20 to -10°C.

Thawing Instructions:

When preparing human serum for use, it should be allowed to thaw gradually, either at room temperature or in a refrigerator at a temperature of 2-8°C overnight to minimize thermodynamic stress to the material.

Appearance:

Straw-like or amber color. It could be slightly cloudy and may contain precipitates due to freezing and/or thawing; this occurrence does not impact culture performance.

Stability:

If stored frozen (-20 to -10°C), product expiration is 5 years from date of manufacture.

Preparation and Use:

We strongly recommend using the entire 100 mL bottle when first thawed. If it is not possible for your intended use, we recommend that you minimize the number of freeze/thaw cycles the product will undergo by aliquoting the serum into containers of a single-use volume.

Biological Testing	
Endotoxin	<10 EU/mL
Hemoglobin	<20.0 mg/dL
Sterility	
Bacteria	Not Detected



Fungi	Not Detected
Mycoplasma	Not Detected
Viral Testing	
HBsAg	Non-Reactive
Anti-HCV	Non-Reactive
Anti-HIV-1/HIV-2	Non-Reactive
Syphilis	Negative
HBV-NAT	Not Detected
HIV-NAT	Not Detected
HCV-NAT	Not Detected
Antibody Screen	Negative
Anti-HBc IgG	Non-Reactive
Anti-HTLV I/II	Non-Reactive
Anti-HBc IgM	Non-Reactive
West Nile Virus NAT	Non-Reactive
Chagas	Non-Reactive
Viral Testing Finished Pool Level	
HBcAB	Non-Reactive
HTLV-I/II	Non-Reactive
West Nile Virus	Not Detected
Chagas	Non-Reactive
Physical Testing	
Osmolality	260-350 mOsm/kg
pH	Test and Report
Biochemistry	
Albumin	Test and Report
ALT (SGPT)	Test and Report
AST (SGOT)	Test and Report
Bilirubin, Total	Test and Report
BUN	Test and Report
Calcium	Test and Report
Chloride	Test and Report
Cholesterol	Test and Report
Creatinine	Test and Report
Glucose	Test and Report
Phosphorus	Test and Report
Potassium	Test and Report
Protein, Total	Test and Report
Sodium	Test and Report
Triglycerides	Test and Report
Uric Acid	Test and Report
Storage	Store frozen at ≤ -10°C
Form	Liquid
Expiration	Five years from the date of manufacture
Handling	Biohazardous; contains human source material. Handle as though capable of transmitting infectious agents.

For Cell Culture, further Manufacturing and Research use only. Not for direct Therapeutic use.