



GOGGLE SANITIZER CABINET

Thomas No. 23A00W341
Mfr. No. GGSN10

Designed to meet the needs of laboratories, the Eisco Goggle Sanitizing Cabinet is safe and easy to use. The cabinet is made from heavy gauge CRC steel and can sanitize up to 36 pairs of goggles in just 15 minutes using Ultraviolet (UV) light.

This process is called Ultraviolet Germicidal Irradiation (UVGI) and has documented safety and efficacy for disinfecting hard surfaces that are directly exposed to UV light.

FEATURES

- Includes several built-in safety features
- Complete UV light sanitation in just 15 minutes
- Accommodates up to 36 pairs of safety goggles or glasses
- Made from heavy gauge steel with a UV resistant enamel finish for durability
- Can be used on a flat surface or hung on the wall (hardware included)
- New and improved packaging for safe and easy shipping



SAFETY FEATURES

- Clearly visible operating instructions on the front of the cabinet (a)
- Viewing window includes a UV filter, which allows users to safely see when the unit is running (b)
- Timer (c) will not begin and UV light will not turn on until doors are CLOSED and LOCKED (d)
- An automatic sensor turns off the UV light if the doors are opened during the sterilization process (e)

HOW IT WORKS

Place up to 36 pairs of safety goggles or glasses into the wire slots in the UV sanitizing cabinet. Once the doors are closed and locked, set the timer for 15 minutes.

High power UV light fills the box and kills all of the germs on the goggles, sterilizing them for the next user.

SPECIFICATIONS

| | |
|----------------------|-------------------------------|
| Input Voltage: | 120V/240V AC +/- 10%, 50/60Hz |
| Input Current: | 0.315A |
| Wattage: | 25 Watts (Max) |
| Size: | 26 x 10 x 29 in |
| Weight: | 54 lbs. |
| Cord Length: | 8 ft |
| Goggle Compartments: | 3.5 x 8 x 3.5 in |