

ENZYCLEAN® II

Dual Enzyme Detergent



Enzyclean® II Dual Enzyme Detergent is specially formulated for use as a presoak or manual detergent. The protease and amylase enzymes in Enzyclean® II break down organic and inorganic soils, residues, films and deposits and ensures thorough and complete cleaning, leaving instruments moving freely and spot-free.

- **Protease enzyme effectively breaks down proteins. Amylase enzyme breaks down carbohydrates.**
- **Reduces likelihood of rust, tarnish or corrosion.**
- **Product safe and compatible with all metals.**

Dual Enzyme Detergent

Order No.	Size	Physical Properties	
B9	1 Gallon (3.8 Liters)	Form:	Liquid
		Color:	Green
		Odor:	Spearmint
		pH Concentration:	8.75

Uses	Temperature	
<input checked="" type="checkbox"/> Presoak	<input type="checkbox"/> Automatic Washer	Usage: 60°-140° F (16°-60° C)
<input checked="" type="checkbox"/> Manual	<input type="checkbox"/> Automatic Endoscope Reprocessor	Storage: Cool, dry area
<input type="checkbox"/> Ultrasonic	<input type="checkbox"/> Cart Washer	

Product Features	Consumer Benefits
Dual Enzyme:	Protease enzyme effectively breaks down proteins. Amylase enzyme breaks down carbohydrates.
Rust Inhibitors:	Reduces likelihood of rust, tarnish or corrosion.
Rinsing Agents:	Rinses spot & residue-free.
Highly Concentrated:	Low cost per use.
Neutral pH:	Product safe and compatible with all metals.
Bacteriostatic:	If stored properly, not favorable to bacterial growth.
Environmentally-Friendly:	Biodegradable, phosphate-free formula.

Procedures	Directions for Use
Directions:	Wear appropriate personal protective equipment based on safe work practices.
Presoak & Manual Cleaning:	1/8 fl. oz. (3.7 mL) to 2 fl. oz. (59 mL) per gallon (3.8 L) of warm water. Allow to soak for 2-3 minutes or longer. Discard solution after each use.
	Dilutions vary depending on temperature, water quality.