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Performance Study of Disposable Sampler (ITM)

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NEST ITM contains guanidine salts and ionic detergents which can lyse viruses and cells to release proteins and nucleic acids. The released proteins, such as DNases and RNases, will be quickly cleaved and denatured, whereas the reducing agent and protective agent in the preservation solution can protect the released nucleic acids from being degraded at normal temperature. ITM can be used for the preservation and inactivation of clinical specimens, such as nasal and pharyngeal swabs. It can effectively prevent changes of microorganisms in such specimens in the types and relative abundance. After being preserved in the virus sampling tube, the nucleic acids in microorganisms can maintain high integrity and can be used for various biological experiments such as enzymatic hydrolysis, PCR and next-generation sequencing. ITM is widely used in the collection and preservation of specimens in hospitals, scientific research institutions, and households.

NEST ITM is designed to be stored at 2-37℃.

Indications for Use

NEST ITM is an enclosed system intended for the collection, stabilization and transportation of pharyngeal and nasal swabs from the collection site to the testing laboratory. the specimen transported in NEST ITM can be used for molecular detection in the laboratory.

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低检出限附近进行测试以研究 ITM 防止核酸降解的能力。The detection limit refers to the lowest concentration of nucleic acids that can be measured in the ITM. In this study, parainfluenza virus type 2, adenovirus and Influenza A were used as research materials and a set of evaluation systems was established through a specific extraction and amplification platform. Tests were conducted under different conditions near the detection limit to study the ability of the ITM to prevent nucleic acids from being degraded.

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Dilute the adenovirus (Guangzhou BDS Biotech Co., Ltd., concentration: 1.57E+5 copies/mL) with HELA cell in a 0.1-fold gradient, add the dilute solution into the ITM at a ratio of 1:10 and then perform extraction with nucleic acid (DNA/RNA) extraction and purification kits (spin column) (SC903-50) (Wuxi Techstar Technology Co., Ltd.). Use the adenovirus nucleic acid test kits (fluorescent isothermal amplification) of Techstar and the Hongshi SLAN-96S device for amplification. Perform the test repeatedly under the concentration of 1.0E + 5copies/mL-1.0E + 4copies/mL for 6 times, and that of 5.0E + 3copies/mL for 20 times. The detection rate should be 100% after 6 repetitions under the concentration of 1.0E+5copies/mL-1.0E+4copies/mL and be at least 95% after 20 repetitions under the concentration of 5.0E+3 copies/mL.

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concentration (copies/mL)	Test1 (Ct)	Test2 (Ct)	Test3 (Ct)	Test4 (Ct)	Test5 (Ct)	Test6 (Ct)	AVERAGE (Ct)	CV
1.0E+5	12.65	12.64	12.85	12.14	12.48	12.95	12.62	2.28%
1.0E+4	15.96	16.81	14.58	15.34	14.96	15.71	15.56	5.07%

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Table 2 Limit test results under the concentration of 5.0E + 3copies/mL of adenovirus

5.0E+3 copies/mL	
Times of repetition	Ct
1	17.58
2	17.28
3	18.05
4	17.62
5	19.34
6	18.09
7	18.39
8	17.39

9	19.25
10	20.39
11	21.85
12	24.39
13	21.36
14	19.27
15	19.57
16	18.46
17	17.69
18	19.2
19	18.08
20	19.36
AVERAGE	19.13
CV	9.25%

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Dilute the Influenza A (Guangzhou BDS Biotech Co., Ltd., concentration: 1.82E+5 copies/mL) with HELA cell in a 0.1-fold gradient, add the dilute solution into the ITM at a ratio of 1:10 and then perform extraction with nucleic acid (DNA/RNA) extraction and purification kits (spin column) (SC903-50) (Wuxi Techstar Technology Co., Ltd.). Use the Influenza A nucleic acid test kits (fluorescent isothermal amplification) of Techstar and the Hongshi SLAN-96S device for amplification. Perform the test repeatedly under the concentration of 1.0E + 5copies/mL-1.0E + 4copies/mL for 6 times, and that of 5.0E + 3copies/mL for 20 times. The detection rate should be 100% after 6 repetitions under the concentration of 1.0E+5copies/mL-1.0E+4copies/mL and be at least 95% after 20 repetitions under the concentration of 5.0E+3 copies/mL.

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Table 4 Limit test results under the concentration of 5.0E + 3copies/mL of Influenza A

5.0E+3 copies/mL		
Times	of	Ct

repetition	
1	23.54
2	30.07
3	26.28
4	29.31
5	30.61
6	27.04
7	26.85
8	25.2
9	28.36
10	25.14
11	29.35
12	26.2
13	25.84
14	24.09
15	22.24
16	25.39
17	26.07
18	30.54
19	28.25
20	32.14
AVERAGE	27.13
CV	9.69%

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Dilute the HPIV-2 virus (Guangzhou BDS Biotech Co., Ltd., concentration: 2.0E+4 copises/mL) with HELA cell in a 0.1-fold gradient, add the dilute solution into the ITM at a ratio of 1:10 and then perform extraction with nucleic acid (DNA/RNA) extraction and purification kits (spin column) (SC903-50) (Wuxi Techstar Technology Co., Ltd.). Use the HPIV-2 virus nucleic acid test kits (fluorescent isothermal amplification) of Techstar and the Hongshi SLAN-96S device for amplification. Perform the test repeatedly under the concentration of 1.0E + 5copies/mL-1.0E + 4copies/mL for 6 times, and that of 5.0E + 3copies/mL for 20 times. The detection rate should be 100% after 6 repetitions under the concentration of 1.0E+5copies/mL-1.0E+4copies/mL and be at least 95% after 20 repetitions under the concentration of 5.0E+3copies/mL.

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concentration (copises/mL)	Test1 (Ct)	Test2 (Ct)	Test3 (Ct)	Test4 (Ct)	Test5 (Ct)	Test6 (Ct)	平均值 (Ct)	CV
1.0E+5	14.39	14.85	14.3	15.2	14.07	14.67	14.58	2.81%
1.0E+4	16.34	16.86	16.53	16.38	16.11	16.04	16.38	1.82%

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Table 6 Limit test results under the concentration of 5.0E + 3copies/mL of HPIV-2 virus

5.0E+3 copises/mL	
Times of repetition	Ct
1	17.45
2	18.4
3	17.6
4	18.1
5	18.23
6	18.15
7	17.94
8	18.04
9	17.36
10	20.5
11	18.64
12	17.93
13	17.37
14	17.62
15	18.05
16	18.54
17	17.33
18	17.11
19	17.95
20	18.12
AVERAGE	18.02
CV	3.9%

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The above test results show that ITM is well consistent when preserving parainfluenza virus type 2, adenovirus and Influenza A under the concentration of 5.0E + 3copies/mL.

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parainfluenza virus type 2, adenovirus and Influenza A were used as research materials and the preservation capacity of ITM was determined by different preservation conditions and duration.

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Dilute the adenovirus to the concentration of $5.0E + 3$ copies/mL with HELA cell, add the dilute solution to the ITM at a ratio of 1:10 and store the solution at the temperature of 4°C and 37°C respectively. On Day 0, 3, 9, 15, 21, and 27 after storage, perform extraction with nucleic acid (DNA/RNA) extraction and purification kits (spin column) (SC903-50) (Wuxi Techstar Technology Co., Ltd.) and use the adenovirus nucleic acid test kits (fluorescent isothermal amplification) of Techstar and the Hongshi SLAN-96S device for amplification. Perform the extraction 20 times repeatedly for each detection point under each storage condition, and the degree of difference in Ct values represents the preservation stability.

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Table 7 Preservation Stability Experiment for adenovirus

Days		0	3	9	15	21	27
4°C	AVERAGE (Ct)	20.24	19.39	22.24	19.25	19.37	20.92
	CV	10.25%	9.28%	10.52%	8.25%	9.3%	12.28%
37°C	AVERAGE (Ct)	21.28	22.92	22.84	19.67	19.6	20.85
	CV	12.2%	10.51%	8.52%	9.71%	9.34%	10.43%

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Dilute the Influenza A to the concentration of $5.0E + 3$ copies/mL with HELA cell, add the dilute solution to the ITM at a ratio of 1:10 and store the solution at the temperature of 4°C and 37°C respectively. On Day 0, 3, 9, 15, 21, and 27 after storage, perform extraction with nucleic acid (DNA/RNA) extraction and purification kits (spin column) (SC903-50) (Wuxi Techstar Technology Co., Ltd.) and use the Influenza A nucleic acid test kits (fluorescent isothermal amplification) of Techstar and the Hongshi SLAN-96S device for amplification. Perform the extraction 20 times repeatedly for each detection point under each storage condition, and the degree of difference in Ct values represents the preservation stability.

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Table 8 Preservation Stability Experiment for Influenza A

Days		0	3	9	15	21	27
4°C	AVERAGE (Ct)	28.49	29.34	31.51	30.51	29.74	30.4

	CV	7.64%	9.25%	10.2%	10.62%	9.71%	9.15%
37°C	AVERAGE (Ct)	30.52	29.84	31.52	31.84	29.54	30.48
	CV	10.25%	12.52%	12.57%	10.23%	9.4%	10.84%

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Dilute the HPIV-2 virus to the concentration of 5.0E + 3copies/mL with HELA cell, add the dilute solution to the ITM at a ratio of 1:10 and store the solution at the temperature of 4 °C and 37 °C respectively. On Day 0, 3, 9, 15, 21, and 27 after storage, perform extraction with nucleic acid (DNA/RNA) extraction and purification kits (spin column) (SC903-50) (Wuxi Techstar Technology Co., Ltd.) and use the HPIV-2 virus nucleic acid test kits (fluorescent isothermal amplification) of Techstar and the Hongshi SLAN-96S device for amplification. Perform the extraction 20 times repeatedly for each detection point under each storage condition, and the degree of difference in Ct values represents the preservation stability.

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Table 9 Preservation Stability Experiment for HPIV-2 virus

Days		0	3	9	15	21	27
4°C	AVERAGE (Ct)	19.04	20.43	19.2	18.83	19.43	20.34
	CV	6.94%	7.54%	6.93%	7.16%	5.39%	7.94%
37°C	AVERAGE (Ct)	22.54	21.56	19.54	20.57	22.63	23.79
	CV	7.93%	6.03%	6.5%	6.15%	5.93%	5.83%

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The above experimental results showed that ITM had good stability in the preservation of adenovirus, influenza A and parainfluenza virus type 2 near the detection limit.

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Three kinds of virus, including adenovirus, Influenza A and parainfluenza virus 2 (10⁷TCID₅₀/ml), were incubated with ITM for 10, 30 and 60 seconds respectively. virus only and ITM only were also incubated accordingly to serve as internal controls. Four days after inoculation, the cells were fixed and stained with 0.06% crystal violet in 1%

glutaraldehyde. Wells that were not stained demonstrated the cytopathic effect (CPE) of the virus. The titer of the virus was recorded as the TCID50.

Inactivation rate:

ITM showed cytotoxicity on MDCK cells (Influenza A), MRC-5 cells (adenovirus) and HELA cells (parainfluenza virus 2) when diluted 1:100 but not cytotoxicity at 1:1,000. The mixture of ITM and virus had similar profiles as ITM alone with 10, 30 and 60 second incubation, ITM rapidly inactivated virus with a >4.0 log reduction in concentration at 10 seconds.

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	10s incubation	30s incubation	60s incubation
Adenovirus only	6.14	6.39	6.45
Adenovirus and ITM	< 2	< 2	< 2
ITM only	< 2	< 2	< 2
	10s incubation	30s incubation	60s incubation
Influenza A only	6.75	6.28	6.69
Influenza A and ITM	< 2	< 2	< 2
ITM only	< 2	< 2	< 2
	10s incubation	30s incubation	60s incubation
Parainfluenza virus 2 only	6.63	6.58	6.1
Parainfluenza virus 2 and ITM	< 2	< 2	< 2
ITM only	< 2	< 2	< 2

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Dilute the adenovirus to the concentration between 1.0E+5 copies/mL-5.0E + 3copies/mL with HELA cell, add the dilute solution to the ITM at a ratio of 1:10. The following commercial nucleic acid extraction reagents were used to determine the compatibility of ITM with nucleic acid extraction reagents.

NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA

- Nucleic acid (DNA/RNA) Extraction Kit (Wuxi Techstar Technology Co., LTD.)
- Nucleic acid (DNA/RNA) Extraction Kit (Da An Gene Co.,Ltd. Of Sun Yat-sen University)
- EasyPure® Viral DNA/ RNA Kit (Beijing TransGen Biotech Co., LTD)
- AxyPrep Body Fluid Viral DNA/RNA Miniprep Kit (Corning Life Sciences(Wujiang)Co. LTD)

The 200µL solution was extracted with the above extraction reagent, and eluent with 100µL elution solution and amplified with the adenovirus nucleic acid detection kit (thermostatic fluorescence method) of Wuxi Techstar Technology Co., LTD. ABI 7500 and Shanghai Hongzhi SLAN-96s were used as detection instruments. Each concentration is repeated 20 times.

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Table 11 Experimental results of extraction platform and amplification instrument compatibility

Concentration of adenovirus	Instrument	Nucleic acid (DNA/RNA) Extraction Kit (Wuxi Techstar Technology Co., LTD.) Positive Results/Total Tests	Nucleic acid (DNA/RNA) Extraction Kit (Da An Gene Co.,Ltd. Of Sun Yat-sen University) Positive Results/Total Tests	EasyPure® Viral DNA/ RNA Kit (Beijing TransGen Biotech Co., LTD) Positive Results/Total Tests	AxyPrep Body Fluid Viral DNA/RNA Miniprep Kit (Corning Life Sciences(Wujiang)Co. LTD) Positive Results/Total Tests
1.0E+5 copies/mL	ABI 7500	20/20	20/20	20/20	20/20
	SLAN-96S	20/20	20/20	20/20	20/20
1.0E+4 copies/mL	ABI 7500	20/20	20/20	20/20	20/20
	SLAN-96S	20/20	20/20	20/20	20/20
5.0E+3 copies/mL	ABI 7500	19/20	19/20	19/20	19/20
	SLAN-96S	20/20	19/20	20/20	19/20

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Dilute the Influenza A to the concentration between 1.0E+5 copies/mL-5.0E + 3copies/mL with HELA cell, add the dilute solution to the ITM at a ratio of 1:10. The following commercial nucleic acid extraction

NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA

reagents were used to determine the compatibility of ITM with nucleic acid extraction reagents.

- Nucleic acid (DNA/RNA) Extraction Kit (Wuxi Techstar Technology Co., LTD.)
- Nucleic acid (DNA/RNA) Extraction Kit (Da An Gene Co., Ltd. Of Sun Yat-sen University)
- EasyPure® Viral DNA/ RNA Kit (Beijing TransGen Biotech Co., LTD)
- AxyPrep Body Fluid Viral DNA/RNA Miniprep Kit (Corning Life Sciences(Wujiang)Co. LTD)

The 200µL solution was extracted with the above extraction reagent, and eluent with 100µL elution solution and amplified with the Influenza A nucleic acid detection kit (thermostatic fluorescence method) of Wuxi Techstar Technology Co., LTD. ABI 7500 and Shanghai Hongzhi SLAN-96s were used as detection instruments. Each concentration is repeated 20 times.

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Table 12 Experimental results of extraction platform and amplification instrument compatibility

Concentration of Influenza A	Instrument	Nucleic acid (DNA/RNA) Extraction Kit (Wuxi Techstar Technology Co., LTD.) Positive Results/Total Tests	Nucleic acid (DNA/RNA) Extraction Kit (Da An Gene Co., Ltd. Of Sun Yat-sen University) Positive Results/Total Tests	EasyPure® Viral DNA/ RNA Kit (Beijing TransGen Biotech Co., LTD) Positive Results/Total Tests	AxyPrep Body Fluid Viral DNA/RNA Miniprep Kit (Corning Life Sciences(Wujiang)Co. LTD) Positive Results/Total Tests
1.0E+5 copies/mL	ABI 7500	20/20	20/20	20/20	20/20
	SLAN-96S	20/20	20/20	20/20	20/20
1.0E+4 copies/mL	ABI 7500	20/20	20/20	20/20	20/20
	SLAN-96S	20/20	20/20	20/20	20/20
5.0E+3 copies/mL	ABI 7500	19/20	19/20	19/20	19/20
	SLAN-96S	19/20	19/20	19/20	19/20

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Dilute the parainfluenza virus 2 to the concentration between 1.0E+5 copies/mL-5.0E + 3copies/mL with HELA cell, add the dilute solution to the ITM at a ratio of 1:10. The following commercial nucleic acid extraction reagents were used to determine the compatibility of ITM with nucleic acid extraction reagents.

- Nucleic acid (DNA/RNA) Extraction Kit (Wuxi Techstar Technology Co., LTD.)
- Nucleic acid (DNA/RNA) Extraction Kit (Da An Gene Co.,Ltd. Of Sun Yat-sen University)
- EasyPure® Viral DNA/ RNA Kit (Beijing TransGen Biotech Co., LTD)
- AxyPrep Body Fluid Viral DNA/RNA Miniprep Kit (Corning Life Sciences(Wujiang)Co. LTD)

The 200µL solution was extracted with the above extraction reagent, and eluent with 100µL elution solution and amplified with the parainfluenza virus 2 nucleic acid detection kit (thermostatic fluorescence method) of Wuxi Techstar Technology Co., LTD. ABI 7500 and Shanghai Hongzhi SLAN-96s were used as detection instruments. Each concentration is repeated 20 times.

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Table 13 Experimental results of extraction platform and amplification instrument compatibility

Concentration of Influenza A	Instrument	Nucleic acid (DNA/RNA) Extraction Kit (Wuxi Techstar Technology Co., LTD.) Positive Results/Total Tests	Nucleic acid (DNA/RNA) Extraction Kit (Da An Gene Co.,Ltd. Of Sun Yat-sen University) Positive Results/Total Tests	EasyPure® Viral DNA/ RNA Kit (Beijing TransGen Biotech Co., LTD) Positive Results/Total Tests	AxyPrep Body Fluid Viral DNA/RNA Miniprep Kit (Corning Life Sciences(Wujiang)Co. LTD) Positive Results/Total Tests
1.0E+5 copies/mL	ABI 7500	20/20	20/20	20/20	20/20
	SLAN-96S	20/20	20/20	20/20	20/20
1.0E+4 copies/mL	ABI 7500	20/20	20/20	20/20	20/20
	SLAN-96S	20/20	20/20	20/20	20/20
5.0E+3 copies/mL	ABI 7500	19/20	20/20	19/20	19/20
	SLAN-96S	19/20	19/20	20/20	20/20

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The above test results indicate that the extraction method and amplification instruments used may not affect the LoD of the downstream assay. The extraction methods listed above all indicate they are compatible with NEST ITM. A similar result is observed when the ABI 7500 and SLAN-96S are compared to each other. But there are many kinds of extracting reagents and testing platforms, Other extraction platform or amplification instrument be used may need additional validation to ensure robust results are provided.

3.2 保存效果对比 Comparison of preservation performance

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The adenovirus with a concentration of 1.0E+4 copies/mL and 5.0E+3 copies/mL diluted with HELA cell was added to ITM at a ratio of 1:10, and the preservation solution was placed at 37℃ for preservation. After 20 and 23 days of storage, the nucleic acid (DNA/RNA) extraction kit of Wuxi Techstar Technology Co., LTD was used for extraction and purification (SC903-50). The amplification reagent and instrument were respectively The adenovirus nucleic acid detection kit (thermostatic fluorescence method) and Hongshi SLAN-96S instrument. Each detection point was repeated 20 times, with the positive detection rate representing the performance of the preservation solution.

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Table 14 Comparison of the test device with predicate device(Adenovirus)

Times	solutions	1.0E+4 copies/mL Positive Results/Total Tests	5.0E+3 copies/mL Positive Results/Total Tests
20 days	ITM	20/20	19/20
	Predicate Device	20/20	20/20
23 days	ITM	20/20	19/20
	Predicate Device	20/20	19/20

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The Influenza A with a concentration of 1.0E+4 copies/mL and 5.0E+3 copies/mL diluted with HELA cell was added to ITM at a ratio of 1:10, and the preservation solution was placed at 37℃ for preservation. After 20 and 23 days of storage, the nucleic acid (DNA/RNA) extraction kit of Wuxi Techstar Technology Co., LTD was used for extraction and purification (SC903-50). The amplification reagent and instrument were respectively The Influenza A nucleic acid detection kit (thermostatic fluorescence method) and Hongshi SLAN-96S instrument. Each detection point was repeated 20 times, with the positive detection rate representing the performance of the preservation solution.

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Table 15 Comparison of the test device with predicate device(Influenza A)

Times	solutions	1.0E+4 copies/mL Positive Results/Total Tests	5.0E+3 copies/mL Positive Results/Total Tests
20 days	ITM	20/20	19/20
	Predicate Device	19/20	20/20
23 days	ITM	19/20	19/20

Predicate Device	19/20	19/20
------------------	-------	-------

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The parainfluenza virus 2 with a concentration of $1.0E+4$ copies/mL and $5.0E+3$ copies/mL diluted with HELA cell was added to ITM at a ratio of 1:10, and the preservation solution was placed at 37°C for preservation. After 20 and 23 days of storage, the nucleic acid (DNA/RNA) extraction kit of Wuxi Techstar Technology Co., LTD was used for extraction and purification (SC903-50). The amplification reagent and instrument were respectively The parainfluenza virus 2 nucleic acid detection kit (thermostatic fluorescence method) and Hongshi SLAN-96S instrument. Each detection point was repeated 20 times, with the positive detection rate representing the performance of the preservation solution.

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Table 16 Comparison of the test device with predicate device(Parainfluenza virus 2)

Times	solutions	$1.0E+4$ copies/mL Positive Results/Total Tests	$5.0E+3$ copies/mL Positive Results/Total Tests
20 days	ITM	20/20	19/20
	Predicate Device	20/20	20/20
23 days	ITM	20/20	20/20
	Predicate Device	20/20	20/20

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The above test results indicate that compared to the predicate device, EST ITM stored RNA at 37°C for 23 days showed no significant difference.