

Technical Data

Bordet Gengou Agar Base

Bordet Gengou Agar Base is recommended for the detection and isolation of *Bordetella pertussis* and *Bordetella parapertussis*. Also used for the "cough plate" method in case of whooping cough.

Composition**

Ingredients	Gms / Litre
Potatoes, infusion from	125.000
Peptic digest of animal tissue	10.000
Sodium chloride	5.500
Agar	20.000
Final pH (at 25°C)	6.7 ± 0.2
**Formula adjusted, standardized to suit performance parameters	

Directions

Suspend 40.00 grams in 1000 ml distilled water containing 10 ml glycerol. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C and aseptically add 15 - 20% sterile, fresh defibrinated blood (sheep, rabbit, human or horse). For selectivity aseptically add rehydrated contents of 2 vials of Bordetella Selective supplement (FD004). Mix thoroughly, taking care to avoid incorporation of air bubbles and pour into sterile Petri plates.

Principle And Interpretation

Bordet Gengou Agar Media were originally formulated by Bordet and Gengou (1) for cultivation of *Bordetella* species. *Bordetella pertussis* is the causative agent of whooping cough and with the help of cough-plate technique, *B. pertussis* can be isolated from pharyngeal extracts, nasopharyngeal secretions and pre-nasal swabs. Kendrick and Eldering (2) modified the original media by replacing 50% human or rabbit blood with 15% sheep blood to make the medium more enriched for detection of *B. pertussis* by the virtue of its haemolytic reaction. Enrichment of the basal media with 25% human blood aids in the detection of Mycobacterium species from small sputum inocula and in Streptomycin sensitivity testing (3).

The medium is highly nutritious thus supports luxuriant growth of *Bordetella* species and can also be used for mass cultivation of *B. pertussis* for vaccine production (2) and for maintaining stock cultures (1).

Potato infusion and peptic digest of animal tissue serve as carbon and nitrogen source while glycerol and blood enrichment provides additional nutrients. Sodium chloride maintains osmotic equilibrium. Incubation should be carried out in a moist chamber (60% humidity) at 37°C for upto 7 days. Medium should not be over dried before use. After 40 hours *B. pertussis* colonies appear smooth, raised, glistening with a zone of haemolysis. Some strains of *Bordetella* are not haemolytic. For confirmation, serodiagnosis and biochemical test should be performed. This medium can be made more selective for *Bordetella*, by using antibiotics like penicillin (4), methicillin (5), cephalexin (6) of which, cephalexin was found to be superior. Cephalexin suppresses unwanted nasopharyngeal growth and significantly increases the isolation rate of *Bordetella* species. Cephalexin is used at a concentration of 40 mg/litre (FD004). Amphotericin B (10 µg/ml) can be added as an antifungal agent to the medium.

For isolation of *B. pertussis* from specimens, use standard procedures. Incubate the plates in a moist chamber at 35-37°C for 7 days and examine daily with or without dissecting microscope (oblique illumination) to detect the presence of *B. pertussis*. Sometimes the accompanying mold colonies can mask the *B. pertussis* colonies. Use sterile scalpel or needle to remove the portion of the agar that contains spreading colonies of moulds. *B. pertussis* colonies may not be visible without the aid of a microscope for 2-4 days. After 7 days of incubation plates may be discarded as negative. Some *Haemophilus* species will grow on *Bordetella* isolation media and cross-react with *B. pertussis* antisera. It may be prudent to rule out X and V factor dependence.

M175

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 2.0% agar gel.

Colour and Clarity of prepared medium

Basal Medium : Light yellow coloured clear to slightly opalescent gel. After addition of glycerol and 15% v/v sterile defibrinated blood: Cherry red coloured opaque gel forms in Petri plates.

Reaction

Reaction of 4% w/v aqueous solution at 25°C. pH : 6.7±0.2

pН

6.50-6.90

Cultural Response

M175: Cultural characteristics observed with added Glycerol and 15% v/v sterile defibrinated blood and Bordetella Selective Supplement (FD004), after an incubation at 35-37°C for 3-4 days.

Organism	Inoculum (CFU)	Growth	Recovery	Haemolysis
Cultural Response				
Bordetella bronchiseptica ATCC 4617	50-100	good-luxuriant	>=50%	gamma
Bordetella parapertussis ATCC 15311	50-100	good-luxuriant	>=50%	gamma
Bordetella pertussis ATCC 8467	50-100	good-luxuriant	>=50%	beta
Staphylococcus aureus ATCC 25923	>=103	inhibited	0%	

Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Reference

- 1. Bordet and Gengou, 1906, Ann. Inst. Pasteur, 20:731.
- 2. Kendrick and Eldering, 1934, Am. J. Public Health, 24:309
- 3. Tarshis M. S. and Frisch A. W., 1951, Am. J. Clin. Pathol., 21:101.
- 4. Flemming A., 1932, J. Path. Bacteriol., 35:831.
- 5. Broome C. V., Fraser D. W. and English J. W., 1979, Internat. Symp. on Pertussis, DHEW, Washington, D.C., 19.
- 6. Suitcliffe E. M. and Abbott J. D., 1972, B. M. J., iii:732.

Revision : 1 / 2011

CE

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516, Swastik Disha Business Park, Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com