



iClean® Viral Transport System (VTM-RT)

Instructions for Use

Intended Use

iClean® Viral Transport System (VTM-RT) is intended for the collection, storage, and transport of clinical specimens containing viruses, chlamydiae, mycoplasma, or ureaplasma from the collection site to the testing laboratory. It provides swabs for sample collection and labeled screw-cap tubes containing Viral Transport Medium (VTM) for culture, storage, and transport of these organisms. iClean® Viral Transport System (VTM-RT) has not been reviewed by the FDA, is intended to serve as a culture media, non-propagating transport, and is to be used by trained and qualified professionals only.

Description and Summary

The laboratory diagnosis of infectious diseases caused by viruses, chlamydiae, mycoplasma, and ureaplasma require proper specimen collection and transport, which can be accomplished with the iClean® Viral Transport System (VTM-RT). iClean® Viral Transport System (VTM-RT) is room temperature stable and consists of a balanced buffer to maintain a neutral pH, antibiotics to inhibit bacterial or fungal contamination, and a source of protein for stabilization.

Hank's balanced salts better preserve the integrity of virus and to improve the extraction from the media. Phenol red is a pH indicator which serves as a visual quality control mechanism. Gentamicin sulfate and amphotericin B inhibit growth of competing bacteria and yeast. The medium is isotonic and non-toxic to mammalian host cells.

Media Formulation

Hanks balanced salt solution(HBSS)
Bovine Serum Albumin (BSA)
Gentamicin sulfate
Amphotericin B
Colistin
L-glutamic acid
HEPES buffer
Phenol red

Storage Conditions & Shelf-Life

iClean® Viral Transport System (VTM-RT) should be stored in a clean, dry, ventilated environment at 2-25°C. The shelf-life of iClean® Viral Transport System (VTM-RT) is 12 months after the manufacture date.

Usage

Proper specimen collection plays a critical role for successful isolation and identification of infectious organisms. Once the sample is collected it should be placed immediately into the media inside the tube and transported to the laboratory. Although the media can maintain even fragile organisms for long periods of time at room temperature, it is recommended that specimens be refrigerated at 2-8°C while in transit. For long term storage, specimens should be frozen at °C-70 or colder. Refer to guidelines and standards for more information on specimen collection procedures.^{2,4,7}

Sterilization

Swabs are sterile. All components are validated and sterilized following ANSI/AAMI/ISO 11137:2006, Sterilization of health care products — Radiation.



Quality Control

Each lot of media is tested for microbial contamination and toxicity to host cells. Refer to publications by ASM^{1,2,4} and CLSI.^{5,6}

Instructions for Use

1. Peel open the sealed pouch containing the swab(s).
2. Remove swab(s) from the pouch and collect a specimen using proper and recommended medical specimen collection procedures for nasopharyngeal and/or oropharyngeal swabs. Do not bend the swab.
3. After the specimen is collected, aseptically remove the cap from the tube and insert the swab into the media.
4. Break the swab shaft at its break-point by bending the swab against the tube wall. After breakage, discard the portion of the shaft outside the tube.
5. Replace the tube cap and secure tightly.
6. Record patient's information on the label.
7. Transport the sample to laboratory for immediate analysis.

Limitations and Precautions

- For In Vitro Diagnostics use only.
- To be used by trained and qualified professionals only.
- This system has not been reviewed by the FDA.
- Intended to serve as a culture media, non-propagating transport.
- If possible, freezing of specimens should be avoided.
- Freezing and thawing may decrease the recovery of viable organisms.
- Improper storage of Sample Collection Kit may lead to decreased antibacterial and antimicrobial activity.
- Do not use wooden shaft swabs with media because they may contain toxins and formaldehydes.
- Using product with tubes or swabs obtained from other sources may affect the performance of the product.
- Do not freeze or heat before use.
- Do not use after expiration date.
- Do not use if product appears damaged, broken, contaminated, is leaking, or the pouch is ripped or open.
- Do not use if the color of the medium has changed to light orange or yellow-pink.
- Do not contact reagent in the tube prior to collecting the sample.
- This is a single-use medical device.
- Do not reuse opened or previously used products.
- Do not ingest the medium.
- Do not attempt to repackage or re-sterilize the product.
- After collection, it is not recommended to remove the swab from the tube.
- Record information on the label prior to delivery and laboratory analysis.

Performance Characteristics

The performance of iClean® Viral Transport System (VTM-RT) was determined by testing of viability of viruses, chlamydiae, mycoplasma, and ureaplasma. 100µL of diluted organism were directly inoculated in triplicate onto swabs. The swabs were then transferred into the transport medium system and held for 0, 24, 48h at both 4°C and RT (20-25°C). At the appropriate time interval, each sample was vortexed and an aliquot of the suspension was inoculated into shell vials or suitable culture media. The cultures were processed by standard laboratory techniques and examined under a specified incubation time. Viability of viruses and chlamydia were determined by fluorescing foci counts and mycoplasma and ureaplasma by CFU counts.

Viruses tested were *Adenovirus*, *Cytomegalovirus*, *Echovirus Type 30*, *Herpes Simplex, Virus Type 1*, *Herpes Simplex, Virus Type 2*, *Influenza A*, *Parainfluenza 3*, *Respiratory Syncytial Virus*, *Varicella Zoster Virus*, *Chlamydia pneumoniae*, *Chlamydia trachomatis*, *Mycoplasma hominis*, *Mycoplasma pneumoniae*, *Ureaplasma urealyticum*.

Results obtained will largely depend on proper and adequate specimen collection, as well as timely transport and processing in the laboratory.

Table 1

Organism	Organism Concentration ^a	Percent Infectivity of Host Cells	Storage Time (hours)	Incubation Time Prior to Reading (hours)	Viability Challenge at Foci of infected cells/200 μ L ^b	
					4°C	RT
Adenovirus	1: 100	2%	0	24	415	415
			24		627	500
			48		795	497
	1: 500	3%	0	24	196	196
			24		229	189
			48		202	64
Cytomegalovirus	1: 10	100%	0	24	822	822
			24		235	50
			48		327	353
	1: 100	100%	0	24	249	249
			24		147	52
			48		121	317
Echovirus Type 30	1: 100	64%	0	24	147	147
			24		515	553
			48		664	799
	1: 500	2.91%	0	24	111	111
			24		328	322
			48		283	192
Herpes Simplex Virus Type 1	1: 10	6%	0	24	285	285
			24		854	432
			48		847	1015
	1: 100	48%	0	24	168	168
			24		127	89

			48		110	142
Herpes Simplex Virus Type 2	1: 10	47%	0	24	139	139
			24		72	110
			48		140	9
	1: 100	97%	0	24	32	32
			24		40	40
			48		23	14

^aThe test swab tip was inoculated in 100µL of suspension and then placed into test tubes with 3mL of iClean transport medium

^bAverage of triplicate tests performed on 200µL of VTM-RT at each time point

Table 2

Organism	Organism Concentration ^a	Percent Infectivity of Host Cells	Storage Time (hours)	Incubation Time Prior to Reading (hours)	Viability Challenge at Foci of infected cells/200µL ^b	
					4°C	RT
Influenza A	1: 50	10%	0	24	384	384
			24		566	339
			48		268	134
	1: 100	12%	0	24	316	316
			24		165	298
			48		118	43
Parainfluenza 3	1: 10	3%	0	48	617	617
			24		40	836
			48		127	163
	1: 100	25%	0	48	445	445
			24		34	352
			48		60	77
Respiratory Syncytial Virus	1: 10	76%	0	24	159	159
			24		196	184
			48		102	157

	1: 100	100%	0	24	31	31
			24		89	67
			48		78	78
Varicella- Zoster Virus	1: 10	100%	0	24	177	177
			24		304	255
			48		46	164
	1: 100	100%	0	24	177	177
			24		109	100
			48		156	143

^aThe test swab tip is inoculated into 100µL of suspension and then placed into test tubes with 3mL of iClean transport medium

^bAverage of triplicate tests performed on 200µL of VTM-RT at each time point

Table 3

Organism	Organism Concentration ^a	Percent Infectivity of Host Cells	Storage Time (hours)	Incubation Time Prior to Reading (hours)	Viability Challenge at Foci of infected cells/200µL ^b	
					4°C	RT
Chlamydia Pneumoniae	1: 10	100%	0	48	202	202
			24		426	524
			48		422	401
	1: 100	100%	0	48	71	71
			24		187	169
			48		134	164
Chlamydia Trachomatis	1: 10	100%	0	48	290	290
			24		283	824
			48		246	336
	1: 100	100%	0	48	83	83
			24		72	188
			48		76	124

^aThe test swab tip is inoculated into 100µL of suspension and then placed into test tubes with 3mL of iClean transport medium

^bAverage of triplicate tests performed on 200µL of VTM-RT at each time point

Table 4

Organism	Organism Concentration ^a	Incubation Time Prior to Reading (Days)	Storage Time (hours)	Viability Challenge at CFU/100µL ^c	
				4°C	RT
Mycoplasma hominis	1: 500	3	0	TNTC ^b	TNTC ^b
			24	TNTC ^b	39
			48	TNTC ^b	86
	1: 1000	5	0	213	213
			24	145	35
			48	179	14
Mycoplasma pneumoniae	Neat	6	0	TNTC ^b	TNTC ^b
			24	TNTC ^b	TNTC ^b
			48	TNTC ^b	1235
	1: 10	6	0	1221	1221
			24	869	337
			48	903	197
Organism	Organism Concentration ^a	Incubation Time Prior to Reading (Days)	Storage Time (hours)	Viability Challenge at CFU/100µL ^c	
				4°C	RT
Ureaplasma urealyticum	1: 500	5	0	TNTC ^b	TNTC ^b
			24	TNTC ^b	TNTC ^b
			48	TNTC ^b	TNTC ^b
	1: 1000	5	0	1122	1122
			24	1379	1081
			48	700	620

^aThe test swab tip is inoculated into 100µL of suspension and then placed into test tubes with 3mL of iClean transport medium

^bTNTC, Too numerous to count

^c Average of triplicate tests performed on 100µL of VTM-RT at each time point

Bibliography

1. Isenberg, H.D., 2004. Clinical Microbiology Procedures Handbook, 2nd. ASM, Washington, DC.
2. Isenberg, H.D., 1998. Essential Procedures for Clinical Microbiology. Chapter 14.12, Page 787. Packaging and shipping infectious Substances. ASM, Washington, D.C.
3. National Committee for Clinical Laboratory Standards (NCCLS). 2003. Quality Control of Microbiological Transport Systems. Approved Standard M40-A, CLSI, Wayne, PA.
4. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.L. Landry and M.A. Pfaller. 2007. Manual of Clinical Microbiology. 9th ed. ASM, Washington, D.C.
5. CLSI. Quality Control of Microbiological Transport Systems; Approved Standard—Second Edition. CLSI document M40-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
6. CLSI. Collection, Transport, Preparation, and Storage of Specimens for Molecular Methods; Approved Guideline. CLSI document MM13-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.
7. 42CFR72. Code of Federal Regulations, Title 42, Volume 7, Part 72. Interstate Shipment of Etiologic Agents.

Manufacturer



Huachenyang (Shenzhen) Technology Co., Ltd.

Area A, F8/Bldg 4, Hengchangrong High Tech Park, Shangnan East Road, Hongtian, Shajing, Bao'an, Shenzhen, 518125, Guangdong CHINA

Website: www.chenyanglobal.com

Tel: 86-755-27393226 Fax: 86-755-27381080

Instruction Approval and Revision Date

Approval Date: 2020-08-18

Revision Date: 2020-08-19

Date of Issue: 2020-08-19

Product Descriptions and Catalog Number

Catalog No.	Description	Pack
CY-F005-20	3ml transport medium vial with iClean® nasal flocked Swab	Single pack, 100 units/box
CY-F005-20	3ml transport medium vial with iClean® buccal flocked Swab	Single pack, 100 units/box
CY-F005-20	3ml transport medium vial with iClean® throat flocked Swab	Single pack, 100 units/box
CY-F005-20	3ml transport medium vial with iClean® nasal flocked Swab and throat flocked Swab	Single pack, 100 units/box
CY-F005-20	3ml transport medium vial with polyester tip swab	Single pack, 100 units/box
CY-F005-20	3ml transport medium vial with 2 polyester tip swabs	Single pack, 100 units/box