

# **Technical Data**

## **Anaerobic Blood Agar Base**

M1345

Anaerobic Blood Agar Base is recommended for isolation and cultivation of Group A and Group B Streptococci from throat cultures and other clinical samples.

### Composition\*\*

Ingredients	Gms / Litre
Casein enzymic hydrolysate	14.500
Papaic digest of soyabean meal	5.000
Sodium chloride	5.000
Growth Factors	1.500
Agar	14.000
Final pH ( at 25°C)	7.3±0.2

<sup>\*\*</sup>Formula adjusted, standardized to suit performance parameters

#### **Directions**

Suspend 40 grams in 990 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure(121°C) for 15 minutes. Cool to 45-50°C. Aseptically add rehydrated contents of 1 vial of Neomycin Supplement (FD149), and 5% v/v sterile defibrinated sheep blood. Mix well and pour into sterile Petri plates.

#### **Principle And Interpretation**

Group B streptococcus (GBS) infection is a common bacterial infection that is rarely serious in adults, but can be life-threatening to newborns. Group A Streptococci commonly causes strep throat and rarely, a potentially deadly destruction of flesh. Anaerobic Blood Agar Base with Neomycin Supplement is used for the isolation of Group A and Group B Streptococci from clinical specimens (1). This medium was originally formulated by Blanchette and Lawrence (2), by addition of the antibiotic Neomycin to sheep blood agar. This addition improved the detection of Group A & B Streptococci, while inhibiting the growth of the other accompanying haemolytic organisms. Casein enzymic hydrolysate and papaic digest of soyabean meal in the medium provide carbon and nitrogenous compounds. Growth factors and defibrinated sheep blood together supply enrichment for growth of fastidious organisms. Sodium chloride helps in maintaining the osmotic equilibrium. Addition of Neomycin supplement (FD149) helps to suppress the normal flora thereby enhancing recovery of Group A and Group B Streptococci.

### **Quality Control**

#### **Appearance**

Cream to yellow homogeneous free flowing powder

#### Gelling

Firm, comparable with 1.4% Agar gel.

#### Colour and Clarity of prepared medium

 $Basal\ medium: Yellow\ coloured\ clear\ to\ slightly\ opalescent\ gel.\ After\ addition\ of\ 5\%v/v\ sterile\ defibrinated\ blood: Cherry\ red\ coloured\ ,\ opaque\ gel\ forms\ in\ Petri\ plates$ 

#### Reaction

Reaction of 4.0% w/v aqueous solution at 25°C. pH: 7.3±0.2

#### pН

7.10-7.50

#### **Cultural Response**

M1345: Cultural characteristics observed in presence of 5-10% CO2 with added 5% v/v sterile defibrinated sheep blood and Neomycin Supplement (FD149), after an incubation at 35-37°C for 24-48 hours.

Organism Inoculum Growth Recovery Haemolysis (CFU)

#### **Cultural Response**

HiMedia Laboratories Technical Data

Escherichia coli ATCC	50-100	none-poor	<=10%	none
Staphylococcus aureus ATCC 25923	50-100	none-poor	<=10%	none
Streptococcus agalactiae ATCC 13813	50-100	good-luxuriant	>=50%	beta
Streptococcus pyogenes	50-100	good-luxuriant	>=50%	beta

#### **Storage and Shelf Life**

Store below  $30^{\circ}\text{C}$  in tightly closed container and the prepared medium at 2 -  $8^{\circ}\text{C}$ . Use before expiry date on label.

#### Reference

Disclaimer:

- 1. Murray P. R., Baron J. H., Pfaller M. A., Jorgensen J. H. and Yolken R. H., (Ed.). 2003, Manual of Clinical Microbiology, 8th Ed., American Society for Microbiology, Washington, D.C.
- 2. Blanchette and Lawrence, 1967, Am. J. Clin. Pathol., 48-411

Revision: 1 / 2011

CE

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.