

Heated Modules

Intelligent, pre-calibrated adapters



A range of intelligent Heated Modules for use with the Vitl Ther-Mix and Flexi-Therm.

Description

Vitl has developed a range of Intelligent Heated Modules for use with the Ther-Mix heated mixer and Flexi-Therm dry block heater.

Each Heated Module is pre-calibrated to ensure accurate temperature control with no need for thermometers or thermocouples and processing extra data.

Modules are available for a range of consumable types such as microtubes, microplates and larger tubes – see www.vitlproducts.com for the full range and specifications.

We can also design bespoke Modules based on the requirements of your method.

Key Features

- Modules can be used on either the Flexi-Therm or Ther-Mix interchangeably
- Pre-calibrated
- Independent control of sample and lid temperatures to reduce condensation and evaporation
- Memory function for previous Module settings
- Direct consumable contact on modules with a lid provides more accurate and efficient sample heating
- Self-limited mixing speed and temperature for sample care and user safety

Related products



Flexi-Therm
V103001



Ther-Mix
V102001



See the whole range on our website or contact sales@vitlproducts.com for more information

T: +44 1233 638 383 **UK**

T: +1 804 381 0905 **USA**