



For collection, transport,  
maintenance and long term freeze  
storage of Viruses, Chlamydia,  
Mycoplasma and Ureaplasma

# SMART TRANSPORT MEDIUM AND ACCUSWAB

Viruses, Chlamydia, Mycoplasma  
and Ureaplasma



## Smart Transport Medium (STM)



1.5 ml Vial  
Phenol Red



1.5 ml Vial  
Clear



3 ml Vial  
Phenol Red



3 ml Vial  
Clear

## AccuSwab™



Nasopharyngeal




Oropharyngeal


## ADVANTAGES

- ▶ Stable at room temperature.
- ▶ Patient sample material and particles are rapidly released and dispersed.
- ▶ Our STM includes antibiotics to help inhibit bacteria and fungal flora in the sample for the most accurate result.
- ▶ The molded break points are both safe and reliable.
- ▶ Our STM is a stand-up and shatter proof tube enabling quality centrifugation of the sample.


## STM15 - 1.5 ml - Phenol Red

	Package	50 units per package
	Vial	Screw cap with tube
	Active ingredient	Phenol Red (11%)
	Chemical name	Phenolsulfonphthalein
	Storage temperature	5 to 25 C
	Intended use	Culture: Viruses, Chlamydia, Micoplasma and Ureaplasma


## STM15 - 1.5 ml - Clear

	Package	50 units per package
	Vial	Screw cap with tube
	Storage temperature	5 to 25 C
	Intended use	Culture: Viruses, Chlamydia, Micoplasma and Ureaplasma


## STM30 - 3 ml - Phenol Red

	Package	50 units per package
	Vial	Screw cap with tube
	Active ingredient	Phenol Red (11%)
	Chemical name	Phenolsulfonphthalein
	Storage temperature	5 to 25 C
	Intended use	Culture: Viruses, Chlamydia, Micoplasma and Ureaplasma


## STM30- 3 ml - Clear

	Package	50 units per package
	Vial	Screw cap with tube
	Storage temperature	5 to 25 C
	Intended use	Culture: Viruses, Chlamydia, Micoplasma and Ureaplasma

## AccuSwab™ Nasopharyngeal

	Package	1000 units per box
	Tip	Nylon Fiber
	Rod	Nylon
	Sterility	Gamma Sterile

## AccuSwab™ Oropharyngeal

	Package	1000 units per box
	Tip	Nylon Fiber
	Rod	Nylon
	Sterility	Gamma Sterile

(AccuSwab and Smart Transport System (STM) are sold separately)