

United Scientific™ Biohazard Bags, PP



These bright red, easy to open, polypropylene (PP) Biohazard Bags are for disposal of biohazard waste and plastic laboratory items such as used tips and tubes. Bags have a clearly printed biohazard symbol and pass ASTM 1922 tear resistance and ASTM 1709 dart impact tests. Steam autoclavable for sterilization of waste before disposal.

- **Polypropylene bags include black biohazard symbol and precautionary procedures printed in four languages (English, Spanish, French, German) for added safety**
- **High temperature steam autoclavable to 135 °C (275 °F)**
- **Pass ASTM 1922 Tear Resistance and ASTM 1709 Dart Impact Tests**
- **55 microns (2.2 mil) thickness**



UNU40106 Biohazard Bag



Easy To Open

Tear-Resistant

Autoclavable

United Scientific™ Biohazard Bags, PP

Item #	Capacity	Dimensions	Weight (approx.)	Quantity per pack
UNU40105	9 Liters	12" W x 18" L	15 g	200
UNU40106	12 Liters	12" W x 24" L	20 g	200
UNU40103	16 Liters	18" W x 22" L	29 g	200
UNU40107	17 Liters	16" W x 24.5" L	24 g	200
UNU40104	37 Liters	25" W x 35" L	64 g	200
UNU40114*	32 Liters	24" W x 31.5" L	46 g	100
UNU40120*	74 Liters	42" W x 42" L	128 g	100

*COMING SOON - Available December 2022



**United
Scientific®**

1001 Technology Way
Libertyville, IL 60048

www.unitedsci.com
847-336-7556