BRANSON

Bransonic® Ultrasonic Cleaning Bath Accessories

General Description

Bransonic accessories are designed to maximize the effectiveness of your Bransonic baths. It is important that parts do not rest on the bottom of the bath, as it can inhibit the ultrasonic activity and possibly damage the unit.

A variety of perforated trays, mesh baskets, beakers, and support racks give you the option of direct or indirect cleaning. All five models come with covers to allow the cleaning bath to heat up more quickly.

Direct Cleaning Method

- Suspends parts using a perforated tray, mesh basket, or a support rack.
- Permits sound waves to penetrate directly to the parts.

Indirect Cleaning Method

- Suspends parts using a solid tray, or beakers with a beaker cover.
- Ideal for cleaning small parts.
- Discarded soil cleans easily from a beaker or solid tray.
- Indirect method is ideal for sample preparation.

Consult with your Branson representative about the Bransonic accessories and solution best suited for your application needs. **Visit us at www.Bransonic.com**



Additional Features

- Perforated trays, solid insert trays, and wire mesh baskets have built-in handles specially designed to prevent accessories from touching the tank bottom.
- Support racks allow total volume of the tank to be used while keeping parts, fixtures, or containers about ½" above the tank bottom.
- Support racks also are useful as countertop drain racks, holding parts above the countertop surface.



BRANSON

Bransonic Accessories – Order using part numbers listed under model numbers

			* 1			
	Ordering Part No.	Model 1800	Model 2800	Model 3800	Model 5800	Model 8800
	Cover (Supplied with Unit)	100-032-515	100-032-517	100-032-519	100-032-510	000-917-050 (SS)
	Perforated Tray	100-410-160	100-410-162	100-410-164	100-410-166	100-410-168
3	Mesh Basket	100-916-333	100-916-334	100-916-335	_	_
and the second	Solid Insert Tray	100-410-170	100-410-172	100-410-174	100-410-176	100-410-178
1881	Beaker Positioning Cover	(1x600 ml) 100-032-516	(2x600 ml) 100-032-518 (2x250 ml) 100-032-527	(2x600 ml) 100-032-520 (3x250 ml) 100-032-528	(4x600 ml) 100-032-511	(6x600 ml) CPN-246-014 (SS)
	Support Racks	CPN-916-039	CPN-916-040	CPN-916-041	CPN-916-042	CPN-916-043
	Beakers (All Models)	250 ml Glass – 000-140-001 400 ml PP – 000-265-061 600 ml Glass – 000-140-004 600 ml SS – 000-410-055				

Americas

Branson Ultrasonics Corp.
41 Eagle Road
Danbury, CT 06810, USA
T: 203-796-0400
F: 203-796-2250
Toll Free: 800-732-9262
www.bransonultrasonics.com

Europe

Branson Ultraschall Niederlassung der Emerson Technologies GmbH & Co. OHG Waldstrasse 53-55 63128 Dietzenbach, Germany T: +49-6074-497-0 F: +49-6074-497-199 www.branson.eu

Δsia

Branson Ultrasonics (Shanghai) Co., Ltd. 758 Rong Le Dong Road Song Jiang, Shanghai, PRC, 201613 T: 86-21-3781-0588 F: 86-21-5774-5100 www.branson.com.cn

