

 BIOFIL®	Document Number	Revision	Page	Document Number
	ATT001	A0	1 of 2	12-03-2020
PRODUCT DATA SHEET		Robotic Tips		



Robotic Tips

Description: Beckman-Sytle Robotic tips, Rack

Purpose: Accurate liquid handling

Materials

Tips: PP (Polypropylene)





Rack box: PP (Polypropylene)

Features

- JET Robotic Tips are made from pure Polypropylene
- Available in both black conductive tips and clear non-conductive tips
- Black conductive tips: liquid level sensing
- Compatibility: Beckman FX/NX®, Multimek AP96® and Biomek 3000® workstations
- Tip format: 96 tips/Rack
- Tip volume range selection: 20µL, 50µL and 250µL
- Available with or without aerosol-resistant filters
- Low Retention/Standard surface
- DNase/RNase-free and non-pyrogenic
- Endotoxin-Safe
- Sterilized by irradiation SAL10⁶ (ISO11137)
- Shelf life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system



Dimensions (Unit: mm)



Robotic Tips

Cat. No.	Max Volume	Tip Color	Surface Type	Filter	Sterilized	Packaging Spec.
ATB000020	20µL	Natural	Standard	No	No	96 tips/rack, 10 racks/pack, 5 packs/case.
AMB000020			Low Retention	No	No	
ATB001020			Standard	No	Yes	
AMB001020			Low Retention	No	Yes	
ATB101020			Standard	Yes	Yes	
AMB101020			Low Retention	Yes	Yes	

 BIOFIL ®	Document Number	Revision	Page	Document Number
	ATT001	A0	2 of 2	12-03-2020
PRODUCT DATA SHEET		Robotic Tips		

Robotic Tips

Cat. No.	Max Volume	Tip Color	Surface Type	Filter	Sterilized	Packaging Spec.
AUB000020	20µL	Black	Standard	No	No	96 tips/rack, 10 racks/pack, 5 packs/case. (Conductive black)
ANB000020			Low Retention	No	No	
AUB001020			Standard	No	Yes	
ANB001020			Low Retention	No	Yes	
AUB101020			Standard	Yes	Yes	
ANB101020			Low Retention	Yes	Yes	
ATB000050	50µL	Natural	Standard	No	No	96 tips/rack, 10 racks/pack, 5 packs/case.
AMB000050			Low Retention	No	No	
ATB001050			Standard	No	Yes	
AMB001050			Low Retention	No	Yes	
ATB101050			Standard	Yes	Yes	
AMB101050			Low Retention	Yes	Yes	
AUB000050		Black	Standard	No	No	96 tips/rack, 10 racks/pack, 5 packs/case. (Conductive black)
ANB000050			Low Retention	No	No	
AUB001050			Standard	No	Yes	
ANB001050			Low Retention	No	Yes	
AUB101050			Standard	Yes	Yes	
ANB101050			Low Retention	Yes	Yes	
ATB000250	250µL	Natural	Standard	No	No	96 tips/rack, 10 racks/pack, 5 packs/case.
AMB000250			Low Retention	No	No	
ATB001250			Standard	No	Yes	
AMB001250			Low Retention	No	Yes	
ATB101180			Standard	Yes	Yes	
AMB101180			Low Retention	Yes	Yes	
AUB000250	250µL	Black	Standard	No	No	96 tips/rack, 10 racks/pack, 5 packs/case. (Conductive black)
ANB000250			Low Retention	No	No	
AUB001250			Standard	No	Yes	
ANB001250			Low Retention	No	Yes	
AUB101180			Standard	Yes	Yes	
ANB101180			Low Retention	Yes	Yes	

Drawn	Approved
Date 03 Dec 2020	Date 03 Dec 2020
Name Alex Fang	Name James