

## Specification Sheet

**PRODUCT:** Hydrochloric Acid, 0.1N

CATALOG NO.: 56095

LOT NO.: TBD

MANUFACTURE DATE TBD

EXPIRATION DATE: 18 Months from the date of manufacture

STORAGE CONDITIONS: 4 - 40°C

### Test Results

Test	Specification	Test Method	Result
Assay	0.099N - 0.101N	Titration	TBD
UV	No Peaks 270-1100nm	UV/Vis Spectrophotometer	TBD
Sterility	No Growth	USP<71>	TBD

Not for direct human, animal or therapeutic use. For Research and Manufacturing Uses.

### Manufacturing Site

Lifecycle Biotechnologies, 5858 Wright Drive, Loveland, Colorado 80538