

StepMate Multi-Dispenser



Optimized for smooth dispensing of even the smallest volumes

Five easy-to-select settings control dispensing volume, right or left-hand operation

By pushing the slide upwards, the syringe is filled up

The stop lever secures the syringe inserted



Main Data Description

Dial setting	1	2	3	4	5		
Number of doses	48	23	15	11	8		
Syringe volume (ml)	Sample volume (µl)					A %	P %
0.50	10	20	30	40	50	± 0.8	≤ 0.7
1.25	25	50	75	100	125	± 0.8	≤ 0.5
2.50	50	100	150	200	250	± 0.8	≤ 0.4
5.00	100	200	300	400	500	± 0.4	≤ 0.3
12.50	250	500	750	1000	1250	± 0.3	≤ 0.3
25.00	500	1000	1500	2000	2500	± 0.2	≤ 0.3
50.00	1000	2000	3000	4000	5000	± 0.2	≤ 0.2

Ordering information

Thomas No.	Cat.No.	Description
1226W05	H721008	StepMate Repeat Pipettor

Features:

- Lightweight and ergonomic
- One-hand operation
- Factory calibrated
- Volumes from 1µl to 5ml can be dispensed
- Maintenance-free
- Equipped with durable tip insertion lever
- Works with eight sizes of disposable, polypropylene syringes from 0.05ml to 50ml



DRAGONLAB