

**Thomas**  
Scientific

**Now Available**  
from Thomas Scientific

# stretchease

Premier Powder-Free Nitrile Examination Gloves

The most comfortable,  
most stretchable glove

**PROPRIETARY  
FORMULA!**

**you've ever worn!**



## SPECIFICATIONS

- Color: periwinkle blue
- Finish: textured fingers
- Type: ambidextrous / non-sterile
- Cuff: beaded
- Length: 9 ½ inches
- Thickness: 3.1 mil (palm), 3.9 mil (fingers), 2.8 mil (cuff),  $\pm 0.01$  mm
- Tensile strength: 14 MPa (min)
- Elasticity/elongation: 500% (min)

## QUALITY STANDARDS

- 1.5 AQL
- Complies with FDA biocompatibility guidance for medical devices
- ISO 13485 certification
- Passes viral penetration test (ASTM F1671)
- 510K/approved for medical use

- Powder free
- More chemical resistant

- Exceptional strength
- Textured



A great alternative for people who are sensitive to latex!

Thomas No.	Mfr. No.	Description	Qty/Box	Qty/Case
1159G03	G4160	Stretchease, Powder Free, Nitrile Exam Goves, XS	100	1,000
1159G04	G4161	Stretchease, Powder Free, Nitrile Exam Goves, S	100	1,000
1159G05	G4162	Stretchease, Powder Free, Nitrile Exam Goves, M	100	1,000
1159G06	G4163	Stretchease, Powder Free, Nitrile Exam Goves, L	100	1,000
1159G07	G4164	Stretchease, Powder Free, Nitrile Exam Goves, XL	100	1,000
1159G08	G4160-250	Stretchease, Powder Free, Nitrile Exam Goves, XS	250	2,500
1159G09	G4161-250	Stretchease, Powder Free, Nitrile Exam Goves, S	250	2,500
1159G10	G4162-250	Stretchease, Powder Free, Nitrile Exam Goves, M	250	2,500
1159G11	G4163-250	Stretchease, Powder Free, Nitrile Exam Goves, L	250	2,500
1159G12	G4164-230	Stretchease, Powder Free, Nitrile Exam Goves, XL	230	2,300

**Larger Case  
Quantity  
Available!**

**Thomas**  
Scientific

ThomasSci.com 833.544.SHIP (7447)  
CustomerService@ThomasSci.com

Connect With Us:

