

Product Name:	ABS VS 5
Description:	<p>ABS VS Class Vapor Shipper</p> <p>Includes 1 Canister</p> <p>Total Vial Capacity (5/Cane): 48</p> <p>Liquid Nitrogen Capacity: 4.3 Liters</p> <p>½ CC Straws (1 Level Bulk): 210</p> <p>½ CC straws (10/Cane): 120</p> <p>Static Evaporation Rate (Liters/day): 0.2</p> <p>Static Holding Time (days): 21</p> <p>Neck Diameter (mm): 51</p> <p>Exterior Height (mm): 492.8</p> <p>Exterior Diameter (mm): 221</p> <p>Canister Height (mm): 279.4</p> <p>Canister Diameter (mm): 46</p> <p>Optional protective shipping container ensure upright shipping</p> <p>High strength neck tube reduces liquid nitrogen loss</p> <p>Lid with locking provision</p> <p>Hydrophobic liquid nitrogen absorbent system</p> <p>2-hour charging capability</p> <p>One-year parts and labor warranty</p> <p>One-year vacuum warranty</p>