

Product Name:	ABS ET 36
Description:	<p>ABS ET Series</p> <p>Includes 6 Canisters</p> <p>Total Vial Capacity (5/Cane): 150</p> <p>½ CC Straw (1 Level Bulk): 780</p> <p>½ CC Straw (10 Cane): 540</p> <p>Static Evaporation Rate (Liters/day): 0.10</p> <p>LN2 Capacity: 36.5</p> <p>Normal Working Time: 224 Days</p> <p>Neck Diameter (mm): 55.4</p> <p>Exterior Height (mm): 690.1</p> <p>Exterior Diameter (mm): 462.3</p> <p>Canister Height (mm): 279.4</p> <p>Canister Diameter (mm): 38.1</p> <p>Designed for small capacity storage</p> <p>Durable aluminum tank with powder coated exterior</p> <p>High strength neck tube reduces liquid nitrogen loss</p> <p>Lid lock provision</p> <p>Low liquid nitrogen consumption</p> <p>Liquid nitrogen measuring rod for level validation</p> <p>One-year warranty on all parts and labor</p> <p>Three-year vacuum warranty</p>