

# safeForce Glass Disposal Boxes



## Safely handle broken glass

*Featuring NEW Carry Handles!*

- Sturdy corrugated cardboard
- 2 mil thick polypropylene liner
- Two models: Bench and Floor

**Keep the lab safe with  
Thomas glass disposal boxes!**

Constructed from rugged corrugated cardboard and supplied with a 2 mil thick puncture-resistant polyethylene liner, the unit safely contains broken glass and spilled liquids. When full, simply pull the safety cap into place and dispose of the entire unit. Striking graphics make the glass disposal box easy to find.

### Box comes in two versions:

Benchtop (8" x 8" x 10") and Floor (12" x 12" x 27") which hold 5.5 and 40 lbs respectively.

Thomas No.	Description	Quantity
20A00T530	safeForce Glass Disposal Box, Floor Model	6/Pack
20A00T531	safeForce Glass Disposal Box, Bench Model	6/Pack