



# Circulating Water Bath Thermostats TW-2.02 – Steel



## Product Description

The stainless steel TW-2.02 convenient and reliable 8.5 liter bath that consistently deliver outstanding results in the lab.

## Features

- Precise and user-friendly temperature control in the bath
- Automatic water level detection system
- Variable water circulation speed
- External water circulation pump to accurately heat other items
- Thorough water mixing to maintain consistent temperature
- Slim, simplified design for quick and easy disassembling and cleaning

## Specifications

Temperature Range in the Bath	Ambient +3 up to 100 °C
Accuracy of Set Temperature Maintenance in the Bath	±0.1 °C
Max Time of Set Temperature Achievement	50 min
Max Bath Capacity	8.5 L
Temperature Stabilization Above the Environmental	3 °C
Water Circulation Speed Variable	2-5 L/min
Voltage and Frequency	100-240 V; 50/60 Hz
Power Consumption	Up to 1.5 kW
Size of the Device (L x W x H)	265 x 325 x 270 mm
Weight	3.8 kg

