



IPTG[®] Isopropyl- β -D-Thiogalactopyranoside

Biotechnology Grade
PURE, Dioxane free
High Efficiency

COMPARE TO OTHER WELL-KNOWN BRANDS!

Package Size	Thomas	EMD Millipore (420322)	Sigma (I6758)
1g	\$26.32	\$59.00	\$57.90
5g	\$63.06	\$149.00	\$255.50
25g	\$210.66	\$419.00	N/A
100g	\$767.22	\$1,059.00	N/A

IPTG Isopropyl- β -D-Thiogalactopyranoside. An inducer of β -galactosidase used to promote the expression of proteins in cells controlled by the lac operator system. For optimum induction, 0.1mM final IPTG concentration in LB (Luria Broth) media is recommended.

CATALOG

C994F28 | 1g, \$26.32

C99429 | 5g, \$63.06

C994F30 | 25g, \$210.66

C994F32 | 100g, \$767.22



Thomas
Scientific