

AASDHPPT Polyclonal Antibody

Catalog number: 11244-1-AP

Size: 31 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

AASDHPPT; AASD-PPT; CGI-

80; DKFZp566E2346; LYS2; LYS5



human brain tissue were subjected to SDS PAGE followed by western blot with 11244-1-AP(AASDHPPT Antibody) at dilution of 1:1000

Applications

Tested applications: ELISA, WB
Species specificity: Human, mouse, rat; other species not tested.
Positive WB detected in: Human brain tissue
Calculated AASDHPPT MW: 36kd
Observed AASDHPPT MW: 36kd
Recommended dilution: WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP =

Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen: Ag1763
GenBank accession number: BC016728
Gene ID (NCBI): 60496
Full name: Amino adipate-semialdehyde dehydrogenase-phosphopantetheinyl transferase

Product information

Purification method: Antigen affinity purification
Storage: PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.