

Lab Specimen Transport Bags

Laboratory Specimen Transport Bags facilitate collection, storage, identification, and transport of laboratory specimens and cultures within the hospital or surgery center, as well as to off-site facilities and laboratories. Bags contain a dual pouch system which separates and protects samples from the accompanying paperwork.

- **NEW! LabGuard® Branding**, the specimen packaging trusted by hospitals and labs for over 25 years
- **NEW! TearZone® EZ-Open Feature** makes specimen bags faster to open than the zipper when a specimen is in processing
- **NEW! Expanded Temperature / Stat Selector Check Boxes** for even greater versatility when using Medegen's Lab Guard specimen bags
- Available in a variety of standard sizes in clear or colored film with easy to use zip closure, reclosable adhesive, or permanent tamper evident adhesive closure systems
- Typical areas of use: Phlebotomy Labs, Chemistry & Hematology, Pathology, Toxicology, Microbiology, Emergency Room, Treatment Areas, Drug Screening, or wherever Routine Clinical Laboratory Specimen Collections are drawn
- **FDA Registered!**



4914



SBL2X69B5



NB4217



4510

Lab Guard® Lab Specimen Transport Bags — Zip Closure

Lab Guard® brand provides convenient, easy-to-use, cost-saving features you just can't get with any other bag in a range of combinations. Choose below from single zip and liquid tight double zip closures, standard or cost saving destroyable biohazard symbol, TearZone® EZ open feature or NO TearZone, paperwork pouch flap or NO flap, STAT+Temp or just Temp selector boxes and more!

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging		Qty/cs
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
SBL2X46B		4x6	Single Zip	Clear/Blk & Org		X	Standard		no flap	DispBag	10/100	1000
	3670	4x6	Single Zip	Clear/Blk & Org	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
	49-91	6x6	Single Zip	Clear/Blk & Red	X		Standard		w/flap	Zip Bag	10/100	1000
SBL2X69B5		6x9	Double Zip	Clear/Blk & Org		X	Destroyable	TearZone®	no flap	DispBag	20/50	1000
SBL2X69B		6x9	Double Zip	Clear/Blk & Org		X	Destroyable	TearZone®	no flap	DispBOX	2/500	1000
	4075	6x9	Single Zip	Clear/Blk & Org	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
	4914	6x9	Single Zip	Clear/Blk & Red	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
	4914.ORG	6x9	Single Zip	Clear/Blk & Org	X		Standard	TearZone®	no flap	Zip Bag	10/100	1000
	D4078	6x10	Single Zip	Clear/Blk & Org	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
	4510	7x9	Slider Zip	Clear/Blk & Org	X		Standard		w/flap	Zip Bag	10/50	500
SBL2X88B		8x8	Single Zip	Clear/Blk & Org		X	Destroyable	TearZone®	no flap	DispBag	20/50	1000
	49-95	8x8	Single Zip	Clear/Blk & Red	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
SBL2X810B5		8x10	Single Zip	Clear/Blk & Org		X	Destroyable	TearZone®	no flap	DispBag	20/50	1000
SBL2X810B		8x10	Single Zip	Clear/Blk & Org		X	Destroyable	TearZone®	no flap	DispBOX	2/500	1000
	4915	8x10	Single Zip	Clear/Blk & Red	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
	4915.ORG	8x10	Single Zip	Clear/Blk & Org	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
SBL2X912B		9x12	Single Zip	Clear/Blk & Org		X	Standard	TearZone®	no flap	DispBag	20/50	1000
	49-99	9x12	Single Zip	Clear/Blk & Red	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
	4925	10x10	Single Zip	Clear/Blk & Org	X		Standard	TearZone®	w/flap	Zip Bag	10/100	1000
SBL2X1118B		11x18	Single Zip	Clear/Blk & Org		X	Standard		no flap	Zip Bag	20/50	1000
SBL2X1215B		12x15	Single Zip	Clear/Blk & Org		X	Standard		no flap	DispBag	40/25	1000
	4916	12x15	Single Zip	Clear/Blk & Org	X		Standard	TearZone®	w/flap	Zip Bag	5/100	500
MGTZ1215		12x15	Single Zip	Clear/Blk & Red		X	Standard		no flap	Zip Bag	5/100	500
	59-94	12x13x4	Single Zip	Clear/Blk & Red	NONE	NONE	Standard		w/flap	Zip Bag	10/100	1000
	NB4217	15x18	Single Zip	Clear/Blk & Org	NONE	NONE	Standard		no doc pch	Zip Bag	5/100	500
ZLAB-TRANSPORT		12.75x12x6	Single Zip	Clear/Blk & Org	NONE	NONE	Standard		no doc pch	Zip Bag	5/100	500
ZLAB1420		14x20	Single Zip	Clear/Blk & Org	NONE	NONE	Standard		no doc pch	Zip Bag	5/100	500
IP2024B3T		20x24	Single Zip	Clear/Blk & Org		X	Standard		w/flap	Zip Bag	1/100	100