

Product Name:	ABT-HC-LP-12
Description:	<p>Premier Glass Door Laboratory Refrigerator  12 Cu. Ft. capacity  One swing glass door, self-closing, right hinged, not reversible  Microprocessor Temperature Controller  Adjustable operating temperature range: 1°-10°C  Digital Temperature Display and Alarm Module  Sample simulated digital temperature display; .1° Resolution Adjustable  temperature display in Celsius or Fahrenheit  Audible and visual high/low temperature alarms  Sensor error alarm  Power failure alarm  Resettable 1-point Min/Max History  Remote alarm contacts  2 probes; 1 in air and 1 in sample bottle  Alarm validation (audible, visual, and remote alarm)  Alarm mute with ring back  Three adjustable shelves (adjustable in ½" increments)  Forced draft air circulation  Keyed door lock  LED interior lights are safety shielded and switch controlled  Swivel casters (3 ½")  probe access port (3/8")  White powder-coated steel exterior and interior  High density urethane foam cabinet insulation  Hydrocarbon, natural refrigerant (R290)  Two year parts and labor warranty, plus an additional three compressor  parts warranty  Exterior dimensions: 25" W x 29 ¾" D x 62 ¼" H without casters; 65 ¾"  H with casters  Cycle defrost  115V, 60 Hz, 3 Amps, 1/5 HP  ETL listed  Shipping weight: 258 lbs.</p>