



WHEATON[®]

Automatic Crimping Tool

High Performance



High Performance Crimping Tool

Wheaton High Performance Crimping Tool is designed with inter-changeable jaws, so that you only need one tool and purchase the size crimper and decappers needed for your laboratory. The jaw sets are available as crimper and decappers which fit standard lab vials seals with aluminum or two-part aluminum/steel caps. This crimping tool has the ability to store separate programs for different caps and cap sizes. The tools are ergonomically designed to reduce strain and arm injury associated with many manual crimping tools.

Product Highlights	
• Fastest and most powerful Wheaton crimping tool	• Strong enough for all steel and magentic caps
• Designed with external power source and cord (no battery)	• Uses interchangeable jaw sets
• Stores up to 9 programs for different caps and seals	• Optional base and mounting kit
• For flip-off caps, available in 13 mm and 20 mm sizes*	

*Not for use with Flip-up caps

Automatic Crimping Tool – High Performance

High Performance Crimping Tool

Cat. No.	Jaw(s) Included	Size (mm)	Qty/Case
W225730-A	—	—	1
W225731-A	Crimper	20	1
W225732-A	Crimper & Decapper	20	1
W225733-A	Flip Cap Crimper & Decapper	20	1

High Performance Jaw Sets

Cat. No.	Jaw Type	Size (mm)	Qty/Case
W225741	Crimper Jaw Set	8	1
W225751	Crimper Jaw Set	11	1
W225752	Decapper Jaw Set	11	1
W225761	Crimper Jaw Set	13	1
W225762	Decapper Jaw Set	13	1
W225763	Flip Cap Crimper Jaw Set	13	1
W225771	Crimper Jaw Set	20	1
W225772	Decapper Jaw Set	20	1
W225773	Flip Cap Crimper Jaw Set	20	1

Automatic Crimping Tool Stand

Cat. No.	Description	Qty/Case
W225701	Crimping Tool Stand	1



For more information on WHEATON products, please visit www.wheaton.com

WHEATON | 1501 North 10th Street | Millville, New Jersey 08332-2038 | USA
800.225.1437 (U.S. & Canada Only) | 856.825.1100 | 856.825.1368 (F)

WSS_056 01/2016