

Specification Sheet

PRODUCT: Sodium Acetate Buffer, 0.05M, pH 4.5

CATALOG NO.: 56444

LOT NO.: TBD

MANUFACTURE DATE TBD

EXPIRATION DATE: 12 Months from the date of manufacture

STORAGE CONDITIONS: 4 - 40°C

Test Results

Test	Specification	Test Method	Result
Acetate Ion	0.05M \pm 0.005	Titration	TBD
pH	4.5 \pm 0.1	Potentiometric	TBD
UV	No Peaks 270-1100nm	UV/Vis Spectrophotometer	TBD
Sterility	No Growth	USP<71>	TBD

Not for direct human, animal or therapeutic use. For Research and Manufacturing Uses.

Manufacturing Site

Lifecycle Biotechnologies, 5858 Wright Drive, Loveland, Colorado 80538