



PRO Scientific Inc.

99 Willenbrock Road, Oxford, Connecticut 06478, USA
Phone:(203)267-4600 Fax:(203)267-4606



Figure 1. Top - Multi-Gen 7, Bottom - Multi-Gen 7XL



Figure 2. Multi-Gen 7 Circular slots on outer shaft



Figure 3. Multi-Gen 7XL Notched "teeth" slots on outer shaft



PRO Scientific Inc.

99 Willenbrock Road, Oxford, Connecticut 06478, USA

Phone:(203)267-4600 Fax:(203)267-4606




Figure 4. Multi-Gen 7XL Packaging, 12 pack only



Figure 5. Multi-Gen 7 Packaging 12 or 24-pack available

Both Multi-Gen 7 and Multi-Gen 7XL Generators

- Easily connects to all PRO motor units with standard Multi-Gen adapter
- Hands-free ejection limits contamination
- Multi-pack design reduces cross contamination risk,
 - o Process each sample with a clean generator
 - o Eliminate cleaning between samples
- 316SS design allows for the utmost in chemical compatibility
- Autoclavable
- Can process in tubes as small as 1.5ml

Multi-Gen 7XL Difference

- Open slotted tooth generators – ideal for tougher samples
- Longer shaft, approximately 115mm – ideal for samples volumes from 1.5/2ml microtube UP TO 15ml conical tubes!!



PRO Scientific Inc.

99 Willenbrock Road, Oxford, Connecticut 06478, USA

Phone:(203)267-4600 Fax:(203)267-4606

PK-01-01200MGS	Standard MULTI-GEN 7 Homogenizing Package; Includes PRO200 homogenizer 115v (p/n 01-01200), toolkit, MULTI-GEN 7 adapter (p/n 07-07200) and 12-pack of MULTI-GEN 7 generators (p/n 02-070MG-12)
PK-01-02200MGS	Standard MULTI-GEN 7 Homogenizing Package; Includes PRO200 homogenizer 220v (p/n 01-02200), toolkit, MULTI-GEN 7 adapter (p/n 07-07200) and 12-pack of MULTI-GEN 7 generators (p/n 02-070MG-12)
PK-01-01200MGD	Deluxe MULTI-GEN 7 Homogenizing Package;Includes PRO200 homogenizer 115v (p/n 01-01200), toolkit, MULTI-GEN 7 adapter (p/n 07-07200), 12-pack of MULTI-GEN 7 generators (p/n 02-070MG-12), PRO200 stand assembly (p/n 80-00200) and MULTI-GEN 7 ejector clamp (p/n 80-00218)
PK-01-02200MGD	Deluxe MULTI-GEN 7 Homogenizing Package;Includes PRO200 homogenizer 220v (p/n 01-02200), toolkit, MULTI-GEN 7 adapter (p/n 07-07200), 12-pack of MULTI-GEN 7 generators (p/n 02-070MG-12), PRO200 stand assembly (p/n 80-00200) and MULTI-GEN 7 ejector clamp (p/n 80-00218)
PK-01-01200MGP	Premium MULTI-GEN 7 Homogenizing Package;Includes PRO200 homogenizer 115v (p/n 01-01200), toolkit, MULTI-GEN 7 adapter (p/n 07-07200), 24-pack of MULTI-GEN 7 generators (p/n 02-070MG-24), PRO200 stand assembly (p/n 80-00200) and MULTI-GEN 7 ejector clamp (p/n 80-00218)
PK-01-02200MGP	Premium MULTI-GEN 7 Homogenizing Package;Includes PRO200 homogenizer 220v (p/n 01-02200), toolkit, MULTI-GEN 7 adapter (p/n 07-07200), 24-pack of MULTI-GEN 7 generators (p/n 02-070MG-24), PRO200 stand assembly (p/n 80-00200) and MULTI-GEN 7 ejector clamp (p/n 80-00218)
PK-01-01200PMGXL	Premium MULTI-GEN 7XL Homogenizing Package;Includes PRO200 homogenizer 115v (p/n 01-01200), toolkit, MULTI-GEN 7XL adapter (p/n 07-07200) and 12-pack of MULTI-GEN 7XL generators (p/n 02-070MGXL-12)
PK-01-02200PMGXL	Premium MULTI-GEN 7XL Homogenizing Package;Includes PRO200 homogenizer 220v (p/n 01-01200), toolkit, MULTI-GEN 7XL adapter (p/n 07-07200) and 12-pack of MULTI-GEN 7XL generators (p/n 02-070MGXL-12)
02-070MG-12	Multi-Gen 7mm 316 stainless steel generators (12/pk)
02-070MG-24	Multi-Gen 7mm 316 stainless steel generators (24/pk)
02-070MG7XL-12	12-pack of MULTI-GEN 7XL generators
02-070MG7XL-12	12-pack of MULTI-GEN 7XL generators
07-07200	Multi-Gen Adapter