

GEB Robotic Tips for Tecan Genesis Freedom®, Freedom Evo®, Miniprep with Liha®.

Product details at-a-glance

GEB Cat#	Max Volume	Tip Color	Low Retention	Sterilized	Filtered (5.1*3.5)mm	Packaging Spec.
AT01MT2-H-N	1000ul	Natural				96 tips/rack, 16racks/case, 1536tips/case.
AT01MT2-H-N-LB	1000ul	Natural	Low Retention			
AT01MT2-H-NS	1000ul	Natural		Sterilized		
AT01MT2-H-NS-LB	1000ul	Natural	Low Retention	Sterilized		
AF01MT2-H-NS	1000ul	Natural		Sterilized	Filtered	
AF01MT2-H-NS-LB	1000ul	Natural	Low Retention	Sterilized	Filtered	
AT01MT2-H-X	1000ul	Black				96 tips/rack, 16racks/case, 1536tips/case. (Conductive black)
AT01MT2-H-X-LB	1000ul	Black	Low Retention			
AT01MT2-H-XS	1000ul	Black		Sterilized		
AT01MT2-H-XS-LB	1000ul	Black	Low Retention	Sterilized		
AF01MT2-H-XS	1000ul	Black		Sterilized	Filtered	
AF01MT2-H-XS-LB	1000ul	Black	Low Retention	Sterilized	Filtered	
Case dimension L*W*H: 405*285*235 mm.			Gross Weight: 4.2kg			

- ✧ 1000μL Liquid Level Sensing Tips
- ✧ 96 Tip/Hanging Tip Rack format
- ✧ Polypropylene
- ✧ Certified DNase/RNase free
- ✧ Nonpyrogenic
- ✧ Presterile/Nonsterile
- ✧ Filtered (Filtered With aerosol-resistant filters) /NonFiltered
- ✧ Endotoxin-Safe meets USP/FDA standards
- ✧ Comply With USP 3127 (Heavy Metal Contamination) standards
- ✧ Compatibility: Tecan Genesis Freedom®, Freedom Evo®, Miniprep with Liha® workstations.

Drawing: (Unmark tolerance Diameter tolerance: $\pm 0.05\text{mm}$, Height tolerance: $\pm 0.25\text{mm}$)

