

## COG8 Polyclonal Antibody

Catalog number: 12661-1-AP


Size: 55 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

COG8; DOR1; FLJ22315



HeLa cells were subjected to SDS PAGE followed by western blot with 12661-1-AP(COG8 Antibody) at dilution of 1:1000

### Applications

|                          |  |
|--------------------------|--|
| Tested applications:     | ELISA, WB                                    |
| Species specificity:     | Human, mouse, rat; other species not tested. |
| Positive WB detected in: | HeLa cells                                   |
| Calculated COG8 MW:      | 612aa, 68kd                                  |
| Observed COG8 MW:        | 68kd   |
| Recommended dilution:    | WB: 1:500-1:5000                             |

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP =

Immunoprecipitation, FC = Flow cytometry

### Immunogen information

|                           |   |
|---------------------------|---|
| Immunogen:                | Ag3341                                  |
| GenBank accession number: | BC017492                                |
| Gene ID (NCBI):           | 84342                                   |
| Full name:                | Component of oligomeric golgi complex 8 |

### Product information

|                      |   |
|----------------------|---|
| Purification method: | Antigen affinity purification                                       |
| Storage:             | PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C. |