

Lovibond® Water Testing

Tintometer® Group



COD Reagents



US EPA / ISO 15705:2002 Compliant

- Suitable for use in other manufacturer's instruments
- Available in 25 and 150 Count Packages
- 3 Testing Ranges
- Mercury-Free Versions Available

www.lovibond.com

COD Reagents



Available for testing 3 ranges, our COD reagents meet US EPA / ISO15705:2002 testing requirements!

Mercury-free versions are also available for non-reporting purposes.

Lovibond® Tintometer® offers COD reagents in packs of 25 or 150 tubes. The Lovibond® line of VARIO chemistries has been formulated for use in our test kits and instruments. They also have the advantage of being suitable for use in other manufacturer's instruments, test kits and other devices.

Ordering Information

Item	Part Number
Low Range COD Reagent (0 - 150 mg/l)	
25 Count	2 42 07 20
150 Count	2 42 07 25
Mid Range COD Reagent (0 - 1,500 mg/l)	
25 Count	2 42 07 21
150 Count	2 42 07 26
High Range COD Reagent (0 - 15,000 mg/l)	
25 Count	2 42 07 22
150 Count	2 42 07 27
Low Range COD Reagent (0 - 150 mg/l)- Mercury Free	
25 Count	2 42 07 10
Mid Range COD Reagent (0 - 1,500 mg/l) - Mercury Free	
25 Count	2 42 07 11
150 Count	2 42 07 16
High Range COD Reagent (0 - 15,000 mg/l) - Mercury Free	
25 Count	2 42 07 12

MD 100 COD Set

The MD 100 COD Set comes complete with everything the user needs for COD analysis.

Simply add sample to the VARIO COD tube test vial, heat it in the RD 125 Reactor, and analyze sample using the MD 100 Single Parameter Colorimeter for COD.

The COD Set is supplied with:

- MD 100 COD Colorimeter
- 16 mm adapter
- 25 tube tests for VARIO Low Range (0-150 mg/L) COD analysis
- 25 tube tests for VARIO Mid-Range (0-1500 mg/L) COD analysis
- RD 125 Reactor
- Tube test stand (24 tube capacity)



Ordering Information

Item	P/N
MD 100 COD Set	27 61 30
Includes: MD 100 COD Colorimeter (with standard accessories), RD 125 Reactor, 1 each of VARIO LR COD Reagent (25 ct.) and VARIO MR COD Reagent (25 ct.) and a Tube Test Rack (24 tube capacity)	

Tintometer Inc.

Tel: +1 941 756 6410
sales@tintometer.us
www.lovibond.us
USA

Tintometer Brazil

Tel: +55 (11) 3230 6410
sales@tintometer.com.br
www.lovibond.com.br
Brazil

Tintometer GmbH.

Tel: +49 (0) 231/94510-0
sales@tintometer.de
www.lovibond.com
Germany

The Tintometer Ltd.

Tel: +44 1980 664800
water.sales@tintometer.com
www.lovibond.com
UK