

alpha Actinin Polyclonal Antibody

Catalog number: 11313-2-AP


Size: 30 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

ACTN1; actinin, alpha 1, ACTN1,
alpha Actinin, Alpha actinin 1, F
actin cross linking protein, Non
muscle alpha actinin 1



C6 cells were subjected to SDS PAGE followed by western blot with 11313-2-AP (α-Actinin Antibody) at dilution of 1:2000

Applications

Tested applications:	ELISA, WB
Cited applications:	WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Cited species:	Rat
Calculated alpha Actinin MW:	103 kDa
Observed alpha Actinin MW:	100-105kd
Positive WB detected in	C6 cells, human spleen tissue, mouse spleen tissue, NIH/3T3 cells
Recommended dilution:	WB: 1:1000-1:10000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag1859
GenBank accession number:	BC015766
Gene ID (NCBI):	87
Full name:	Actinin, alpha 1

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.