

Exclusive Promotional Offer:

Buy 10 Corning® 50 mL FBS
bottles and get 5 at no additional
cost.

For purchases April 1, 2024 through June 30, 2024



Thomas No.	Corning Cat. No.	Description
CHM00A380	35-010-CF	Corning® Fetal Bovine Serum, 50 mL, Regular, USDA Approved Origin
CHM00A381	35-011-CF	Corning® Fetal Bovine Serum, 50 mL, Regular, USDA Approved Origin (Heat Inactivated)
CHM00A382	35-015-CF	Corning® Fetal Bovine Serum, 50 mL, Premium, United States Origin



- USDA origin is collected at government approved harvest facilities in Mexico and Central America, while US origin is collected at government approved abattoirs within the United States.
- True pooled and triple 0.1 micron sterile filtered that passes sterility, mycoplasma, and 9CFR viral testing requirements.
- United States Origin FBS granted a Certificate of Suitability, R0-CEP 2018-271, by the European Directorate for the Quality of Medicines (EDQM).
- For HI - Validated process used to heat inactivate serum at 56°C ± 2°C for 30 minutes via bulk inactivation of the entire lot in the pooling tank or individual bottles in a water bath.

Scan below to learn more:


To redeem this promotion:

Contact your local Corning Account Manager, with code "lot 50" in the subject line.

Name

Email

Phone



Terms and Conditions: Offer valid in the U.S. and Canada only and void where prohibited by government, law, or company policy. Offer valid for end-user customers only for purchases April 1, 2024 through June 30, 2024 or while supply lasts. Redemption information must be received by July 15, 2024. Please allow 4-6 weeks for delivery.
Reservation requests for qualifying 50mL lots are prohibited.

For more specific information on claims, visit www.corning.com/certificates.
Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not for use in humans.

