

Nasopharyngeal Test Swab (100mm Break Point)**Thomas Number: 21A00F582**

Handle Length (mm) [in]	(152) [5.98]
Handle Diameter (mm) [in]	(2.4) [0.095]
Flock Head Length (mm) [in]	(18.42) [0.73]
Breakpoint (mm) [in]	(100) [3.94]
Individual Sterile Package Dimensions (mm) [in]	(7 x 190) [0.28 x 7.48]
Handle Material	Polypropylene
Flock Material	Polyamide (Nylon)
Sterilization Method	Ethylene Oxide Gas
Approved Jurisdictions	USA, Canada, and EU
Swab Quantity per Shelf Pack	100
Shelf Pack Quantity per Shipper Box	10
Shipper Boxes per Pallet	42
Shipper Box Weight (kg) [lbs]	(3.12) [6.88]
Shipper Box Outer Dimensions (mm) [in]	(505 x 400 x 241) [19-7/8 x 15-3/4 x 9-1/2]