

501017
Lid Assembly with Knob

WARING
COMMERCIAL


SS115

Dry Blending Container

Container Information

500 ml Working Capacity

Unique Blade for Dry
Homogenizing and Pulverizing



501681
Container / stainless steel with Blending Assy.

501682
Blending Assy.

Illustration#	Part#	Description
1	003504	Cap Nut
2	024225	Blade
3	003609	Washer / stainless steel
4	003508	Washer / fairprene
5	010012	Bearing Holder
6	009724	Washer / nylon
7	003537	Washer / stainless steel 3 Req.'d
8	003509	Washer / neoprene (rubber) 4 Req.'d
9	003510	Washer / fiber
10	002406	Socket Head Drive Shaft
11	009780	Hex Nut
12	012860	32 Oz. Stainless Steel Container
13	009762	Gasket
14	015055	Container Adapter