

Slim T™ Mop



Micronova's **Slim T™ Mop** is lightweight and easily slips on the T-Mop™ series of adapters. Slim T covers are excellent for cleaning ceilings and walls, lyophilizer shelving and around equipment where a low profile is required. The Slim T Mop heads also work well on soft wall enclosures.

Slim T covers come in a variety of fabrics and can be used for residue and particle removal, dispersing disinfectants, or other cleaning solutions. Slim T covers can be autoclaved and are available sterile for aseptic environments

- Slim T Mops work well on ceilings, walls, floors, windows, soft wall enclosures, equipment shelves and surfaces.
- The slip-on covers fit on QP, QD and QDST series T-Mop adapters, with 180-degree swivel
- Slim T covers also work with Micronova's Window Washing Tool available in 12" & 18" widths
- Covers available in a range of fabrics with a thin or thick polyester urethane foam interior to match your cleaning requirements
- Available in 10", 12", 14", or 18" lengths
- Compatible with alcohol, strong disinfectants and steam sterilization
- Slim T covers are suitable for ISO 5 to ISO 9 environments.

Manufactured for use in Controlled Environments

- ISO 5 thru 9 Cleanrooms
- ISO 9001 Certified Manufacturing Facility
- Autoclavable (if applicable)
- Good chemical resistance to disinfectants
- Minimum Irradiation Dose: 25 kGy
10⁻⁶ Sterility Assurance Level
- In compliance with AAMI Guidelines for Sterilization (if applicable)



Microfiber and NovaPoly Slim T Mops on 18, 14 and 10 in T-Mop Adapters

NovaPoly, 10 in, Slim T Mop on T-Mop Adapter



PolySorb, 18 in Slim T Mop on T-Mop Adapter



Microfiber, thick, 12 in Slim T Mop on Window Washing Tool

SPECIFICATIONS

Materials of Construction:

Urethane Foam Interior

Cleaning Surfaces:

- NovaPoly™:** 100% Smooth polyester knit provides excellent disinfectant dispersal
- PolySorb™:** 100% Textured polyester assists in removing residues and dispersing disinfectant
- Microfiber:** 80% Polyester, 20% polyamide mop is highly absorbent and effective in dry and wet mopping

Dimensions:

- Width: 2.5 in (6 cm)
- Length: 10 in to 18 in (25 cm to 45 cm)



Micronova Manufacturing, Inc.
3431 W. Lomita Boulevard
Torrance, CA 90505-5010 USA

Tel: 310/784-6990
Fax: 310/784-6980

SalesInfo@micronova-mfg.com
www.micronova-mfg.com