

Product Description

Powered by V-Spin – the world's most advanced, patented mixing and liquid-handling technology – the Intelli-Mixer RM-2 is truly the latest in cutting-edge science. V-spin harnesses the energy of an unequaled stepper motor, allowing for a wide variety of mixing, vibrating, and rotating cycles. This multipurpose unit is a rotator, vortex mixer, and blotter all in one. Designed to be cold room safe, and with an unprecedented 34 modes of operation, this versatile piece is unrivaled in today's market.

Features

- · Rocker, Vortex Mixer, and Blotter
- 34 automated intelligent mixing programs
- 3 Year Warranty

Specifications

Rotation Speed	1 - 99 RPM
Rotation Speed Settings Step	1 RPM
Overall Number of Programs	28 Standard + 3 Custom
Environmental Temperature	10 - 45 °C
Relative Air Humidity (at 20°C)	up to 80 %
Power Supply Adapter	DC 12V, 1.5 A

Accessories



























Order: (747) 226-3330

info@elminorthamerica.com 5115 Douglas Fir Rd Suite G | Calabasas, Caifornia 91302