

## KIMTECH™ Vista™ Nitrile Exam Glove

Dark blue, powder-free, non-sterile, chlorinated, ambidextrous exam glove with textured fingertips and beaded cuff. Packed in dispenser box.

## Product Information

Code	Size	Palm Width* (mm ± 5)	Length* (mm)	Quantity	
				Per Dispenser	Per Case
62825	XS	75	240	200	2000
62826	S	85	240	200	2000
62827	M	95	240	200	2000
62828	L	105	240	200	2000
62829	XL	115	240	180	1800

\*As measured by ASTM D6319, AQL 4.0

Physical Property	Unit	Target	Test Method
Residual powder	mg/glove	≤2.0	ASTM D 6124, AQL 4.0
Watertightness	N/A	Pass	ASTM D 5151, AQL 1.5

Thickness	Target (mm)	Target (mil)	Test Method
Finger	0.11	4.33	ATSM D 6319, AQL 4.0
Palm	0.06	2.36	
Cuff	0.05	1.97	

Property	Unit	Target	Test Method (2.5 AQL) <sup>1</sup>
Tensile strength before aging	Mpa	≥14.0	ATSM D 412
Tensile strength after aging	Mpa	≥14.0	ATSM D 412, ATSM D 573
Ultimate elongation before aging	%	≥500	ATSM D 412
Ultimate elongation after aging	%	≥400	ATSM D 412, ATSM D 573



**Workmanship Requirement:** Gloves shall be made in accordance with reasonable industry practice with respect to defects, dirt and contamination. Gloves should be packaged as specified having described size and color in the appropriate quantity as listed. Shipping Case and inner packaging shall be labeled with product name, size, lot/batch and catalog number.

**For Additional information:** Contact your Kimberly-Clark representative or our Customer Service Team

- North America (888)-346-GOKC (4652), [KCPIInfo@kcc.com](mailto:KCPIInfo@kcc.com)
- Asia-Pacific [kimtech.asia@kcc.com](mailto:kimtech.asia@kcc.com)

**Confidentiality.** The information contained in this document is intended for the use of the addressee(s) only and may contain privileged, confidential, or proprietary information that may not be distributed publicly without Kimberly-Clark's written consent.

**Certificate of Conformance** available online at [www.kimtech.com/certificates](http://www.kimtech.com/certificates).