



separopore®

bioaffinity matrix for protein purification

Separopore® 4B is a high gel strength bioagarose matrix for gel filtration and bioaffinity chromatography. It is ideal for both analytical and preparative protein purifications, providing good separation over a wide molecular weight range.

CATALOG #: **C995Z56** | Separopore® (Agarose) 4B, 100ml, **\$128.00**

CATALOG #: **C995Z57** | Separopore® (Agarose) 4B, 500ml, **\$475.21**

CATALOG #: **C995Z60** | Separopore® (Agarose) 4B-CL (cross-linked), 100ml, **\$170.68**

CATALOG #: **C995Z61** | Separopore® (Agarose) 4B-CL (cross-linked), 500ml, **\$546.90**

