

Specification Sheet

PRODUCT: Potassium Phosphate Buffer, 0.05M, pH 7.20

CATALOG NO.:	56436	LOT NO.:	TBD
MANUFACTURE DATE	TBD	EXPIRATION DATE:	12 Months from the date of manufacture
STORAGE CONDITIONS:	4 - 40°C		

Test Results

Test	Specification	Test Method	Result
Phosphate Ion	0.05M \pm 0.005	Titration	TBD
pH	7.20 \pm 0.05	Potentiometric	TBD
UV	No Peaks 270-1100nm	UV/Vis Spectrophotometer	TBD
Sterility	No Growth	USP<71>	TBD

Not for direct human, animal or therapeutic use. For Research and Manufacturing Uses.

Manufacturing Site

Lifecycle Biotechnologies, 5858 Wright Drive, Loveland, Colorado 80538