

Specialty Bags

UV Protection Zip Bags (Amber Color)

Take comfort in knowing the integrity of critical light sensitive products won't be compromised during testing, shipping and storage. Contents can be viewed through the film. Film contains UV blockers to prevent bag contents from harmful light impact. Premium thick, white zipper opens quickly and easily. Complies with USP procedure 66-1 with light transmission not exceeding 10% at any wavelength. No Print.

Minigrip Code#	Bag Color	Size (inches)	Qty/cs	Packaging
MGUV3P0203	Amber	2 x 3	1,000	100
MGUV3P02509	Amber	2.5 x 9	1,000	100
MGUV3P0305	Amber	3 x 5	1,000	100
MGUV3P0406	Amber	4 x 6	1,000	100
MGUV3P0508	Amber	5 x 8	1,000	100
MGUV3P0608	Amber	6 x 8	1,000	100
MGUV3P0808	Amber	8 x 8	1,000	100
MGUV3P0814	Amber	8 x 14	500	100
MGUV3P1212	Amber	12 x 12	500	100



MGUV3P0814

UV Protection Adhesive Closure Bag (White/Gray Opaque)

Speci-Gard® style Press and Close® Adhesive Closure Bag features opaque gray white film that prevents light from entering the bag, protecting light sensitive contents. Access to contents is available with EZ-Tear Perforation. See page 11, item SPECI-UV.

Anti-Static Zip Bags (Pink Tint Color)

- Bags qualify as Type 2 under MIL-B-81705D entitled barrier material flexible and electrostatic free
- Bags are tested to dissipate an induced 5000 volt charge to technical zero in accordance with FTS 101C, method 4046.1
- Smaller sizes packaged in DispenserBag® convenience packaging
- Pink color film indicates anti-static properties; thick, durable white zip closure. No Print.

Minigrip Code#	Bag Color	Size (inches)	Qty/cs	Packaging
MGST4P020503	Pink Tint	2.5 x 3	1,000	100
MGST4P0305	Pink Tint	3 x 5	1,000	100
MGST4P0406	Pink Tint	4 x 6	1,000	100
MGST4P0608	Pink Tint	6 x 8	1,000	100
MGST4P0609	Pink Tint	6 x 9	1,000	100
MGST4P0810	Pink Tint	8 x 10	1,000	100
MGST4P0912	Pink Tint	9 x 12	500	50
MGST4P1212	Pink Tint	12 x 12	500	50
MGST4P1318	Pink Tint	13 x 18	500	50



MGST4P1212

(continued on pg 18)