

RM637 Lecithin soya, 30%

Product Number Packing

RM637	: 100G
RM637	: 500G

Product Information

Product Code	: RM637
Product Name	: Lecithin soya, 30%
Synonym	: L- α -Lecithin, From Soyabean; L- α -phosphatidylcholine, 30%; Soya lecithin, 30%
CAS No.	: 8002-43-5
EC No.	: 232-307-2
HS Code	: 2923 20 00
Shelf Life	: 4 years

Technical Specification

Appearance	: White to yellow to orange to brown powder
Solubility	: 100 mg soluble in 1 mL of chloroform
Acid value	: ≤ 36.0 mg KOH/g
Heavy metals (as Pb)	: $\leq 0.002\%$
Water (K.F.)	: $\leq 1.00\%$
Phosphatides (as acetone insoluble %)	: $\geq 95.00\%$
Phosphatidylcholine (HPLC)	: min. 30.00 %
Total phosphatidylcholine (HPLC)	: min. 32.00 %

Risk and Safety Information

WGK	: 1
RTECS	: OG7565000
Storage Temperature($^{\circ}$ C)	: Store at 2 - 8 $^{\circ}$ C

Transport Information

Marine Pollutant	: No
ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods