

## PRODUCT DATA SHEET

Standing: 2021-07-21

LAUDA Varioshake VS 8 BE

Shaker 230 V; 50/60 Hz

Part Number: L003056

### Features

- Shaker with bidirectional motion, heavy-duty, low-noise in a compact design
- Electronic speed control, stepless, gentle start-up
- Clear control panel for easy operation with rotary knobs
- Timer adjustable up to 60 min. or continuous operation
- Shaking table made of anodized aluminum with four plastic pins
- Drive by an overload-protected AC motor
- Compact, low-wear mechanism with mass balancing for a secure stand
- Constant speed, independent of load
- Excellent for continuous use
- Outer housing made of electrolytically galvanized sheet steel, powder-coated



### Technical Features

|                            |                                      |
|----------------------------|--------------------------------------|
| Ambient temperature range  | 10 ... 60 °C                         |
| Shaking tray (WxD)         | 330 x 330 mm                         |
| Power consumption max.     | 0.07 kW                              |
| Shaking amplitude          | 20 mm                                |
| Shaking frequency          | 20 ... 300 1/min                     |
| Load max.                  | 8 kg                                 |
| Overall dimensions (WxDxH) | 350 x 375 x 160 mm                   |
| Power plug                 | Power cord with angled plug (CEE7/7) |
| Power supply               | 230 V; 50/60 Hz                      |

Reserve technical changes

LAUDA DR. R. WOBSE GMBH & CO. KG  
Pfarrstraße 41/43 • 97922 Lauda-Königshofen  
Postfach 1251 • 97912 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0 • F + 49 (0) 9343 503-222  
info@lauda.de • www.lauda.de  
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen  
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:  
LAUDA DR. R. WOBSE Verwaltungs-GmbH  
Sitz Lauda-Königshofen  
Registergericht Mannheim • HRB 560226

Geschäftsführer:  
Dr. Gunther Wobse (Vors.), Dr. Mario Englert,  
Dr. Marc Stricker  
Beirat: Dr. Gerhard Wobse